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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Orthosis with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

MAINTENANCE

Do not bleach No chemical cleaning

Do not iron Do not tumble-dry

Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Lumbalgia
 - Sciatic lumbalgia and lumbar cruralgia
 - Spondylarthrosis and lumbar discopathies
 - Slight traumas to the lumbosacral column and paravertebral lumbar contractures
 - Post-operative recovery
 - After-effects of stable compression fractures of lumbar vertebrae
 - Treatment of fracture of transverse processes
- Only REF.1322-1322SH-1354-1321-1353:
- Degenerative lumbar scoliosis
 - Spondylolysis with or without vertebral listhesis

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Base fabric:
 - 98% reinforced anaelastic cotton for optimal containment
 - Breathable and natural for best comfort
 - Exclusive design
- Anatomical design shaped for the legs. Available in two sizes:
 - Regular fit (REF.1355-1354-1353)
 - Tailored fit (REF.1323-1322-1321)
 - Tailored fit, lowered on the front (REF.1322SH)
- Traction system:
 - Reinforced elastic straps for targeted pressure and perfect stabilization of the lumbar rachis, latex free
 - Pulley system for easy fastening
 - Ferrules in a soft polymeric material to facilitate grip and removal of the brace
- Front Velcro® fastening system, with comfortable grip for easy closing. Shaped for maximum comfort.
- 2 rigid front and back thermoformable pads, on medical prescription.
- Optimal rigid stays, moldable on medical prescription, to be inserted in the 4 linear back pockets.
- Adapts to the pathology of the patient.

PUTTING ON THE APPLIANCE

REF.1322-1322SH-1323-1354-1355-1325-1356

- Open up the orthosis and provisionally affix the Velcro® straps to the strips on the side (fig. A). Position the support in such a way that the inside label faces upwards.
- Apply the support and use the middle Velcro® closure to fix in place securely (fig. B). Ensure that the closure is precisely centred.
- Slowly tighten up the support by securing the lower Velcro® straps to the front strips (fig. C).
- Then secure the upper straps to the front Velcro® strips (fig. D).

REF.1321-1353

- Open up the orthosis and provisionally affix the Velcro® straps to the strips on the side. Put on the support in such a way that the central Velcro® strip is positioned precisely between the breast bone and the pubic bone (fig. E).
- Guide the side hip strap on the right of the orthosis over the designated Velcro® closure. To do this, insert fingers into the accompanying pockets (upper loop) and generate the appropriate tension (fig. F).
- Complete securing the side hip strap (lower loop) (fig. G).
- Subsequently tighten the orthosis and secure the lower tensioning straps, one over the other (left over right, or vice versa) at the height of the pubic bone (fig. H).
- Secure the upper tensioning straps in the same way (left over right, or vice versa) and select the required height:
 - Overlapping with the lower tensioning straps (fig. I).
 - On top of the lower tensioning straps.

To apply the rear rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing, fig.L)
- apply the adhesive Velcro® to the sides of the central rib (fig.M)
- attach the pad with Velcro® to the internal dorsal part of the orthosis (fig.N)

To apply the front rigid pad (heat-mouldable):

- apply the adhesive padding to the plastic part (concave side-without any writing)
- attach the pad with Velcro® to the front closure of the orthosis

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