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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Epicondylitis elbow strap with pressure pad

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF.24022
Size	universal
Colour	grey
fits right and left	

MAINTENANCE

- ⌘ Do not bleach
- ⌘ No chemical cleaning
- ⌘ Do not iron
- ⌘ Do not tumble-dry
- ⌘ Washing instructions: Remove the pads. Hand wash the bandage, the cover of the cooling pad and the pressure pad at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Acute or chronic epicondylitis and epitrochleitis of the elbow
- Muscle contractures
- Prevention of inflammatory processes during sports activity

CONTRAINDICATIONS

There are no known contraindications at this time.

FEATURES AND MATERIALS

- Polyamide structure for optimal pressure distribution
- Pressure pad to focus the thrust; made of non-slip TPE for an easy positioning of the brace
- Lycra padding coupled with polyurethane foam
- Adjustable strap for accurate compression control
- Velcro® fastener with retaining ring
- Adjustable strap length with practical Y-shaped Velcro®.

APPLICATION

- 1 Apply the pressure pad in the specific Velcro® area, with the long side perpendicular to the long side of the strap.
- 2 Put the strap on:
 - Flatten it out without pulling it out of the ring
 - Slip it on from your hand and pull it up along your arm (fig.A)
 - Position it 3 or 4 fingers below the epicondyle (fig.B) and rotate it until the pressure pad is aligned with the extensor muscles of the hand (fig.C).
- 3 Tighten the Velcro® strap by pulling until the white line appears inside the ring (fig.D). Tighten without causing discomfort or a tingling sensation.

NOTE: for epitrochleitis, wear the brace by positioning the pressure pad diametrically opposite (fig.E)

ONLY FOR THE FIRST APPLICATION

- 1 If necessary, pull out the strap to adjust the length:
 - Remove the Y-shaped Velcro® (fig.F)
 - Cut off the excess part (fig.G)
 - Re-attach the Y-shaped Velcro® (fig.H)

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