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**PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

## Ankle brace with steel stays and figure-8 elastic straps

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF.24081			
Size	0	1	2	3
Circumference ankle cm	16/19	19/22	22/26	26/30
Colour	grey			

fits right and left

### MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap.
- Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Conservative treatment following a grade 2 sprain trauma
- Acute tendinopathies
- Post-surgical rehabilitation
- Aid to resuming sporting activity after grade 2 and 3 sprain trauma
- Chronic joint pain caused by tibiotarsal arthrosis

### CONTRAINDICATIONS

Currently no known

### FEATURES AND MATERIALS

- Three-dimensional mesh fabric for tightness and maximum breathability
- Additional lateral and medial reinforcement fabric to protect against contact with the stays
- Moldable medial and lateral stabilising stays made of steel rods
- Open heel with Velcro® closure for easy donning
- Drawstring adjustment system with quick-release buckle for fast and precise adaptation of the brace to the morphology of the ankle.
  - Lycra pocket for storing the lanyard
- Elastic crossover straps with Velcro® closure on the instep for maximum stability
- Anatomical and low-profile design to fit most shoes

### PUTTING ON THE APPLIANCE

- 1 Loosen lanyard and tie rods and open the back Velcro®.
- 2 Slide the forefoot into the ankle strap from the back (fig. A).
- 3 Close the ankle strap with the rear Velcro®, taking care to align the latter with the corresponding orange line (fig. B).
- 4 Tighten the ankle strap with the special lanyard (fig.C1) and fasten it with the buckle (fig.C2). Insert the excess lanyard into the Lycra pocket on the front flap (fig.C3).
- 5 Position the medial elastic strap:
  - grasp the medially sewn strap and stretch it over the back of the foot towards the opposite side (fig.D1);
  - loop it around the ankle and Velcro® it above the medial malleolus (fig.D2).
- 6 Position the elastic side strap:
  - grasp the laterally sewn strap and stretch it above the back of the foot towards the opposite side (fig.E1);
  - loop it around the ankle and Velcro® it above the lateral malleolus (fig.E2);
- 7 After removing the brace, remember to widen the lanyard before putting it on again (fig.F1, F2).

### ONLY FOR THE FIRST APPLICATION

- 1 If necessary remove the steel stays(fig.G) and mold them according to the patient's morphology.
- 2 The stays can be inserted in two different positions to avoid compression of the malleoli (fig.H1, H2).

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