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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Post-operative knee orthosis with graduated movement

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. The device should not be worn in the vicinity of open flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use. Contact an orthopaedic technician for the replacement of used parts.

SELECTION/SIZE

Item	REF.24045
Size	universal
Length of Hinges cm	56
Color	grey

Fits right and left

MAINTENANCE

-  Do not bleach  No chemical cleaning
-  Do not iron  Do not tumble-dry
-  Washing instructions: Rods: wash with a sponge soaked in warm water and mild soap; dry with a cloth. Padding: remove the paddings and wash them in warm water with mild soap; allow to dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-operative treatment that requires control of knee flexion-extension ligament reconstruction, anterior tibial apophysis transposition, fractures of the distal third of the femur, fractures of the proximal third of the tibia, fractures of the patella, tendon damage to the extensor apparatus
- Post trauma phase of knee sprains
- Varus or valgus osteotomies

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Lightweight aluminium alloy adjustable hinge with integrated strap loops.
- New „PUSH&SET“ flexion-extension adjustment joint, made of a lightweight polymer material that can be set intuitively and safely:
 - Flexion: from -10° to +120° with intervals of 10°
 - Extension: from -10° to +90° with intervals of 10°
- Rapid locking of the joint with a special button, from -10° to +90° with intervals of 10°
- Removable and washable condyle protection pads
- Pads in MTP to limit migration of guide and protect leg from direct contact with the rods; lined only externally, cuttable, for a good fit around the leg
- Possibility of shortening the orthosis using guide marks engraved on the lightweight

ADAPTING

- Open all straps and temporarily close them upon themselves.
- Separate the padding from the rods with the straps.
- Apply the padding such that it fits closely to the leg, closing it at the back with the appropriate Velcro® fasteners: first the thigh padding, then the calf padding (fig.A). In order to position the padding at the optimal height, place the orthosis next to the leg. If necessary the padding may be shortened to adapt it to the patient's legs (shorten the padding along the thinner parts in such a way that the padding coloses properly).
- Model the rods (by bending them) to fit them to the shape of the patients' legs (fig.B).
- Adjust the flexion extension as described (points 10-12, hinge adjustment):perform the same operation on both joints.
- Fix the rods to the pads (fig.C): the rods should follow the median line of the leg and the joints should be aligned with the patella so as to rest on the condyles. To align the rods with the median line of the leg, use the lateral malleolus and the large trochanter as reference points.
- Ensure that the Velcro® straps are perfectly aligned with the corresponding paddings.
- Ensure that the two joints are at the same height.
- Velcro® -fasten the straps after inserting them in the appropriate rings (fig.D-E): fasten first the thigh straps being careful with the one closest to the knee; then proceed by fastening the calf straps. Tighten the straps so as to avoid vertical migration of the guide.

NOTE: The bending procedure should be performed by hand, with the help of the grips or laying the orthosis on a rounded body (chair, edge of a bed); do not use your own knee as a lever.

REGULATION OF THE JOINTS

- Regulate flexion (fig.F1):
 - Pull the ratchet outwards
 - Turn it until the central notch coincides with the desired degree of deflection
 - Release the pawl, engaging the teeth in the slot
- Regulate extension (fig.F2):
 - Pull the ratchet outwards
 - Turn it until the central notch coincides with the desired degree of extension
 - Release the pawl, engaging the teeth in the slot
- To block the joint (fig.F3):
 - Turn the hinge until the lock button coincides with the desired degree of flexion and extension (orange graduated scale)
 - Use snap fastener to lock button into place, causing the teeth to engage in the slot provided.

RODS SHORTENING PROCEDURE

- If rod length adjustment is required (refer to guide marks indicated in fig.G), proceed as follows:
- Detach the proximal (thigh) or distal (calf) padding from the rods.
 - Disengage the fastening straps from the terminal loops.
 - Position the rod on the edge of a stable work surface, ensuring that the pre-defined cut mark is precisely aligned with the edge. Apply firm downward pressure until the rod fractures cleanly at the mark.
 - Smooth the cut edge using abrasive paper or a file to eliminate any sharp edges or burrs.
 - Re-thread the straps through the terminal loops.
 - Trim and reattach the padding to match the new rod length, ensuring proper anatomical alignment and comfort.

RECOMMENDATIONS FOR THE PATIENT

To easily wash the immobilised limb without removing the brace, it is recommended to:

- Open the upper part (straps and padding) and wash the thigh
- Close the upper part
- Open the lower part (straps and padding) and wash the leg
- Close the lower part
- Tighten the strap to prevent vertical migration of the brace.
- To wash the limb we recommend using a damp sponge.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.