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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Eplistik Plus / Eplistik32: Elasticated fabric elbow support
Manulastik Plus / Manulastik19: Elastic wrist kneep brace
Genulastik Plus / Genulastik04: Mesh elastic knee brace
Malleoelastik Plus / Malleoelastik37: Elastic knit ankle support

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure on the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or (REF-90104F) strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never adjust the adjustment made by the doctor/physiotherapist or orthopedic technician. The orthosis is made and to be used by one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

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MAINTENANCE
☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: Hand wash at max.30°C with neutral soap. Do not dry away from heat sources.)
Do not dispose of the product or any of its components into the environment.

INDICATIONS
Eplistik Plus/Eplistik32 - Lateral epicondylitis (tennis elbow) and medial epicondylitis (golfer's elbow) • Chronic pain caused by bursitis or tendinitis. **Manulastik Plus/Manulastik19** - Wrist and carpal sprains • Tendinitis wrist • Arthrosis and arthritis of wrist

• After wrist fractures • Post-surgery wrist treatment. **Genulastik Plus/Genulastik04** - Knee pain due to bursitis • Early stage arthrosis • Chondropathy/Chondromalacia patella • Mild patella instability • Prevention of patella laxation during sports.

Malleoelastik Plus/Malleoelastik37 - 1st degree tibiotarsal sprain • Peroneal or extensor tendinopathy.

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS
• Stretch knit fabric featuring areas of different targeted compression • Anatomically designed

Eplistik Plus/Eplistik32 - Silicone inserts shaped to redistribute compression on the epicondyle and epitrochlea and pressure pads on the forearm • The edging is strengthened with silicone trim so the shoe fits snugly and keeps its shape when being pulled on.

Manulastik Plus/Manulastik19 - Thumb hole • Polymeric palm strip • Velcro® fastening system • Foam-reinforced for softness, traction resistant during application • Silicone insert to fit the knipr brace process.

Genulastik Plus/Genulastik04 - Silicone patella stabilizer • Soft lateral stays to keep brace correctly positioned (anti-tumbler) • Remoiced cotton siseo in edge, so hanes and traction-resistant during application.

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Eplistik Plus/Eplistik32 1 Slide on the elbow support with the label facing upwards: the larger part of the silicone insert should rest on the forearm. 2 Position the elbow support so that the openings in the silicone insert are centred around the epicondyle and epitrochlea.

Malleoelastik Plus/Malleoelastik37 1 Slide on the wrist support with the label facing upwards. 2 Position the wrist support so that the thumb goes through the hole provided and the stiff strip rests against the palm of the hand.

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.