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EN – The descriptions and images in this document are for example and commercial purposes only. Orthoservice reserves the right to modify them according to its needs.

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OSUK Ltd
Units 478/50 North East BIC
Warefield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200 · Fax 0039 0362 18 58 065
info@roplusten.com · www.roplusten.com

elbo PRO

ORTHOSERVICE
RO+TEN
Take care feel better

ORTHOSERVICE
RO+TEN
Take care feel better

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CUSTOMIZABLE
FIT

EASY
ADJUSTING

 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Telescopic elbow brace with graduated movement

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. In addition, we advise on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds. Contact an orthopaedic technician for the replacement of used parts.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.24017R - REF.24017L
Size	universal
Length of the hinges cm	36-46

Colour: dark gray · Right or left

MAINTENANCE

-  Do not bleach
-  Do not iron
-  Do not tumble-dry
-  No chemical cleaning
-  Do not tumble-dry
-  Washing instructions: **Hinges:** wash with a sponge soaked in lukewarm water and neutral soap; dry with a cloth
- Padding:** remove the padding and wash in lukewarm water with mild soap; leave to dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Compound fractures of elbow, distal humerus, proximal radius and ulna
- Immobilisation post-operative following tendon and ligament repairs
- Useful in the rehabilitation phase following surgical and conservative treatment of elbow fractures
- Elbow subluxations and dislocations
- Chronic elbow injuries
- Elbow hyperextension
- Post-traumatic treatment

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Lateral articulated hinge made of a light aluminium alloy. Telescopic length adjustment via easy-to-use buttons on the valves.
- Polyamide valve with mouldable aluminium arm grip for perfect individual fit.
- New „Pull&Set“ flexion-extension adjustment joint, made of a polymer lightweight material for an intuitive and safe regulation:
 - Flexion: from -10° to +120° with intervals of 10°
 - Extension: -10° to +90° with intervals of 10°
 - fast lock joint with a specific button, from -10° to +90° with intervals of 10°
- Length-adjustable straps
- Strap clips supplied
- Heat-sealed padding made of bi-coupled MTP, with soft lycra in contact with the skin
- Length-adjustable shoulder strap, with padded shoulder protector
- Neutral hand support included

APPLICATION FIRST APPLICATION

- 1 Open all straps and close them provisionally on themselves.
- 2 Place the brace on the limb by centering the joint on the elbow joint (fig. A)
- 3 Regulate the length of the telescopic hinges of the arm:
 - Press the button and slide the external valve to the desired length (fig.B1)
 - Release the button by inserting it in the corresponding hole on the metal part (fig.B2)
 - You can move the internal valve (by two positions) making it click with the aid of a screwdriver (fig.B3)
- 4 Regulate the length of the telescopic hinge of the forearm:
 - Press the button and slide the external valve to the desired length (fig.C1)
 - Release the button by inserting it in the corresponding hole on the metal part (fig.C2)
 - You can move the internal valve (by two positions) making it click with the aid of a screwdriver (fig.C3)
- 5 Shape the aluminium arm sockets by hand according to the circumference of the limb (fig. D).
- 6 Fasten the brace to the arm (fig. E) by passing the straps through the appropriate loops and fastening them with Velcro® (first the distal one, then the proximal one). If necessary, shorten the straps by cutting off the excess and reapplying the Y shaped Velcro® tip.
- 7 Attach the brace to the forearm (fig. F) by passing the straps through the strap loops and fastening them with Velcro® (first the distal one, then the proximal one). If necessary, shorten the straps by cutting off the excess and reapplying the Y shaped Velcro® tip.
- 8 If necessary, the supplied strap fasteners can be applied to improve the self-positioning of the straps themselves: slip them over the aluminium arm sockets (fig.G1) and Velcro® them to the padding (fig.G2).
- 9 Make sure the joint is always centered on the elbow joint
- 10 Regulate flexion (fig.H1):
 - Pull the deflection slider outwards
 - Rotate it until the central notch on the slider coincides with the desired degree of deflection
- 11 Regulate extension (fig.H2):
 - Pull the extension slider outwards
 - Rotate it until the central notch on the slider coincides with the desired degree of extension
 - Release the slider, engaging the teeth in the slot
- 12 To block joint (fig.H3):
 - Rotate the hinge until the fast lock button doesn't coincide with the desired degree of flexion-extension (graduated orange scale)
- 13 Snap lock the button making the teeth engage in the appropriate slot
- 14 If necessary, the supplied shoulder strap with shoulder protector can be attached:
 - Insert the loop in the wrist strap (fig.I1).
 - Carry the shoulder strap over the ipsilateral shoulder (Fig.I2), pass it behind the back, under the contralateral armpit and attach it to the same ring by means of the snap fastener (fig.I3).

NEUTRAL HAND SUPPORT APPLICATION

- 1 Insert the supplied sockets into the slots on the hand support (fig. L1)
- 2 Place the hand support on the first forearm valley to fit sleeves in the corresponding holes (fig.L2)
- 3 Insert the screws in the sleeves and tighten all the way (fig.L3)

PATIENT ADVICE

To easily wash the immobilized limb without removing the brace we advise you:

- 1 To open the upper part (straps and padding) and wash the arm
- 2 Close the upper part
- 3 Open lower part (straps and padding) and wash the forearm
- 4 Close the lower part. Tighten straps to avoid the brace from moving upwards

To wash the limb we suggest using a moistened sponge

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