LÖZDZ, UZMANGI IZLASIET UN SAGLABĀJIET INSTRUKCIJAS

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Telescopic elbow brace with graduated movement

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully, do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds. Contact an orthopaedic technician for the replacement of used parts.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe reaction, contact the manufacturer immediately. Do not use the device in the proximity of free flames or strong electromagnetic fields. Do not use the product if the patient has known intolerances or allergies to one or more components in the product. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF.24017R - REF.24017L
Size	universale
Length of the hinges cm	36-46

Colour: dark gray - Right or left

MAINTENANCE

- ☒ Do not bleach ☒ No chemical cleaning
 - ☒ Do not iron ☒ Do not tumble-dry
 - ☒ Washing instructions: **Hinges:** wash with a sponge soaked in lukewarm water and neutral soap; dry with a cloth
 - Padding:** remove the padding and wash in lukewarm water with mild soap/leave to dry away from heat sources.
- Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Compound fractures of elbow, distal humerus, proximal radius and ulna
- Post-operative immobilisation following tendon and ligament repairs
- Useful in the rehabilitation phase following surgical and conservative treatment of elbow fractures
- Elbow subluxations and luxations
- Chronic elbow injuries
- Elbow hyperextension
- Post-traumatic treatment

CONTRAINDICATIONS

Contraindications known at the time of writing

CHARACTERISTICS AND MATERIALS

- Lateral articulated hinge made of a light aluminium alloy. Telescopic length adjustment via easy-to-use buttons on the valves.
- Polyamide valve with mouldable aluminium arm grip for perfect individual fit.
- New „PulIiSet“ flexion-extension adjustment joint, made of a polymer lightweight material for an intuitive and safe regulation:
 - Flexion: from -10° to +120° with intervals of 10°
 - Extension: -10° to +90° with intervals of 10°
 - fast lock joint with a specific button, from -10° to +90° with intervals of 10°
- Length-adjustable straps
- Strap clips supplied
- Heat-sealed padding made of bi-coupled MTP with soft lycra in contact with the skin
- Length-adjustable shoulder strap, with padded shoulder protector
- Neutral hand support included

APPLICATION FIRST APPLICATION

- 1 Open all straps and close them provisionally on themselves.
- 2 Place the brace on the limb by centering the joint on the elbow joint (fig. A)
- 3 Adjust the length of the telescopic hinges of the arm:
 - Press the button and slide the external valve to the desired length (fig.B1)
 - Release the button by inserting it in the corresponding hole on the metal part (fig.B2)
 - You can move the internal valve (by two positions) making it click with the aid of a screwdriver (fig.B3)
- 4 Regulate the length of the telescopic hinge of the forearm:
 - Press the button and slide the external valve to the desired length (fig.C1)
 - Release the button by inserting it in the corresponding hole on the metal part (fig.C2)
 - You can move the internal valve (by two positions) making it click with the aid of a screwdriver (fig.C3)
- 5 Shape the aluminium arm sockets by hand according to the circumference of the limb (fig. D).
- 6 Fasten the brace to the arm (fig. E) by passing the straps through the appropriate loops and fastening them with Velcro® (first the distal one, then the proximal one). If necessary, shorten the straps by cutting off the excess and reapplying the Y shaped Velcro® tip.
- 7 Attach the brace to the forearm (fig. F) by passing the straps through the strap loops and fastening them with Velcro® (first the distal one, then the proximal one). If necessary, shorten the straps by cutting off the excess and reapplying the Y shaped Velcro® tip.
- 8 If necessary, the supplied strap fasteners can be applied to improve the self-positioning of the straps themselves: slip them over the aluminium arm sockets (fig.G1) and Velcro® them to the padding (fig.G2).
- 9 Make sure the joint is always centered on the elbow joint
- 10 Regulate flexion (fig.H1):
 - Pull the deflection slider outwards
 - Rotate it until the central notch on the slider coincides with the desired degree of deflection
 - Release the slider, engaging the teeth in the slot
- 11 Regulate extension (fig.H2):
 - Pull the extension slider outwards
 - Rotate it until the central notch on the slider coincides with the desired degree of extension
 - Release the slider, engaging the teeth in the slot
- 12 To lock joint (fig.H3):
 - Rotate the hinge until the fast lock button doesn't coincide with the desired degree of flexion-extension (graduated orange scale)
 - Release the button, making it click with the aid of a screwdriver
- 13 Snap lock the button making the teeth engage in the appropriate slot
- 14 If necessary, the supplied shoulder strap with shoulder protector can be attached:
 - Insert the loop in the wrist strap (fig.I1).
 - Carry the shoulder strap over the ipsilateral shoulder (fig.I2), pass it behind the back, under the contralateral arm and fasten it to the wrist strap (fig.I3).

NEUTRAL HAND SUPPORT APPLICATION

- 1 Insert the supplied sockets into the slots on the hand support (fig. L1)
- 2 Place the hand support on the first forearm valve to fit sleeves in the corresponding holes (fig.L2)
- 3 Insert the screws in the sleeves and tighten all the way (fig.L3)

PATIENT ADVICE

- 1 To easily wash the immobilized limb without removing the brace we advise you:
 - 1 To open the upper part (straps and padding) and wash the arm
 - 2 Close the upper part
 - 3 Open lower part (straps and padding) and wash the forearm
 - 4 Close the lower part. Tighten straps to avoid the brace from moving upwards
- 2 Wash the limb we suggest using a moistened sponge

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