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UK PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

### Gonarthritis orthosis

**DECLARATION OF CONFORMITY**  
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**TRADEMARKS OF MATERIALS**  
Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**  
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

**WARNINGS**  
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use. Contact an orthopaedic technician for the replacement of used parts

**SELECTION/SIZE**

Code	REF.24014R - REF.24014L						
Size	XS	S	M	L	XL	XXL	3XL
Circumf. 15cm above knee center cm	33/39	39/45	45/52	52/58	58/66	67/74	74/81
Height cm	40	40	40	40	40	40	40

**Right**= for medial gonarthrosis of the right knee (varus right knee) or for lateral gonarthrosis of the left knee (valgus left knee)  
**Left**= for medial gonarthrosis of the left knee (varus left knee) or for lateral gonarthrosis of the right knee (valgus right knee)

**Colour:** dark gray - Right or left

- MAINTENANCE**
- Do not bleach  No chemical cleaning
  - Do not iron  Do not tumble-dry
  - Washing instructions:
    - Take the stay out
    - Wash in lukewarm water with neutral soap
    - Leave to dry away from heat sources
    - Reinsert the stay

Do not dispose of the product or any of its components into the environment.

- INDICATIONS**
- Moderate unicompartmental gonarthrosis
  - Post-operative cartilage or meniscus surgery

**CONTRAINDICATIONS**

Currently no known

- CHARACTERISTICS AND MATERIALS**
- Patented fastening and tensioning system:
    - interlocking buckles and tensioning levers make it fast and easy
    - Ring unloads brace pressure for better comfort
  - Efficient:
    - Three-point pressure system lightens the load on the painful side of the knee
    - Support when needed: tensioned while walking and in stand-by when seated
  - Breathable elastic fabric:
    - Front opening for an easy fit
    - Soft pocket tips to facilitate fastening
  - Technopolymer valve for optimal pressure distribution and effective anchorage to the limb.
    - Single-sided steel stay with flexion (20°, 40°) and extension (0°, 10°) limitation
    - Polymer stay with integrated unloader ring, for correct positioning of the brace
    - Gel pad for additional protection against pressure from the thrust ring
    - Padding for additional valve protection
    - Velcro® strap adjustment

- APPLICATION FIRST APPLICATION**
- Make sure that the two front straps (orange and white) are positioned in the corresponding service hooks around the kneecap and that the levers are open. **CAUTION:** Do not open the Velcro straps.
  - Open the Velcro® fasteners on the fabric base.
  - Wrap the fabric around the knee, keeping the orange label at the top and Velcro® it closed, using the pockets (fig.A1). Ensure that the opening is centred on the kneecap (fig.A2).
  - Hook the interlocking buckles into the corresponding orange and white slots on the ring (fig.B1,B2)
  - Adjust the lengths of the rear (fig.C1) and front (fig.C2) straps by means of the Velcro® straps so that they are neither too loose nor too tight. If too long, cut off the excess and reapply the Y-shaped Velcro® ends(fig.C3).
  - Tension the straps, closing the lower (fig.D1) and upper (fig.D2) levers. Make sure they snap in.
  - Stand up, take a few steps, check that the pressure on the condyle is not excessive but sufficient to unload the contralateral part of the knee to ease the pain. Proceed with any adjustment by repeating and following the steps above in the same order
  - It is possible to shape the stay based on the patient's needs.
  - To regulate flexion-extension:
    - Remove screws and plastic cover (fig.E1)
    - Insert wedges at the desired degrees (fig.E2)
    - Re-attach plastic cover and screws (fig.E1)

- SUBSEQUENT APPLICATIONS**
- CAUTION** - never open the Velcro® unless you want to readjust the orthosis.
- Make sure that the two front straps (orange and white) are positioned in the corresponding service hooks around the kneecap and that the levers are open.
  - Open the Velcro® fasteners of the fabric base.
  - Wrap the fabric around the knee, keeping the orange label at the top and Velcro® it closed, using the pockets provided (fig.F1). Ensure that the opening is centred on the kneecap (fig.F2).
  - Hook the interlocking buckles into the corresponding orange and white slots on the ring (fig.G1,G2)
  - Tension the straps, closing the lower (fig.H1) and upper (fig.H2) levers. Make sure they snap in.

- ORTHOSIS REMOVAL**
- Open the lever (fig.I)
  - Unlock the buckle and fit it in the corresponding orange (fig. L1) and white ( fig.L2) hook slot around the knee cap .
  - Open the Velcro® fasteners from the fabric base and remove the brace.

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