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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

**Knee brace with polycentric hinged rods for F-E control and accessories for compression and cold therapy**

**DECLARATION OF CONFORMITY**

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

**TRADEMARKS OF MATERIALS**

Velcro® is a registered trademark of Velcro Industries B.V.

**SAFETY PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of fire flames or strong electromagneticfields. Do not apply in direct contact with open wounds.

**WARNINGS**

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**COMPRESSION THERAPY**

DO NOT overinflate or overtighten the brace. Ensure compression does not impair blood flow or fluid circulation.

**COLD THERAPY**

The device may reach temperatures low enough to cause serious injury, including tissue necrosis. Always observe the following precautions:

- Frequently check skin condition under the gel pack. If any skin changes are noticed in the treated area, discontinue use immediately and consult a physician
- NEVER apply the gel pack directly to the skin. Always use a barrier (e.g., included mesh sleeve)
- Do not exceed 30 minutes of continuous use

Do not heat the gel pack (microwave, oven, hot water, etc.)

Do not sit or lie on the gel pack

Do not sleep with the gel pack applied

Limit device use in patients with cognitive impairments, sensory disorders, or any condition affecting pain perception.

Patients under sedatives or medications that impair sensory perception must be closely monitored during use.

Always discontinue use at the first sign of cold-related injury.

<b>SELECTION/SIZE</b>	
<b>Code</b>	<b>REF.26100</b>
<b>Size</b>	<b>Universal</b>
<b>Lenght cm</b>	<b>34</b>
<b>Color</b>	<b>black</b>
<b>Fits right and left</b>	

<b>MAINTENANCE</b>	
<span><span></span><span> </span></span> Do not bleach <span><span></span><span> </span></span> No chemical cleaning	
<span><span></span><span> </span></span> Do not iron <span><span></span><span> </span></span> Do not tumble-dry	
<span><span></span><span> </span></span> Knee brace washing instructions:	
• Hand wash in lukewarm water with mild soap;	
• Air dry away from heat sources	
Gel pack washing instructions:	
• Hand wash in lukewarm water with mild soap;	
• Air dry away from heat sources	

Do not dispose of the device or any of its components in the environment.

**INDICATIONS**

Indicated in all cases requiring physiological stabilization of the knee joint, with optional cooling and compression (to relieve inflammation, swelling, edema, and pain), such as:

- Mild to moderate collateral ligament instability
- Meniscus injuries
- Post-surgical recovery
- Osteoarthritis
- Patellofemoral pain

**CONTRAINDICATIONS**

DO NOT use this product on patients affected by: vascular disorders, diabetes, sickle cell anemia, local circulatory disorders, tissue inflammation, cold hypersensitivity or history of cold-related injuries

**CHARACTERISTICS AND MATERIALS**

• Made of soft dual-lined foam with integrated fastening straps:

- Velcro®-compatible outer layer
- Inner layer in elastic fabric
- Wrap-around design with practical Velcro® closure

• Universal size, thanks to two extension straps for larger limb circumferences, which are included

• Polycentric hinge system easily adjustable externally via pins:

- Flexion settings 0°,10°, 20°, 40°, 60°, 90°, 150°

- Extension settings 0°,10°, 20°, 40°, 60°, 90°

- Included accessories for cold and compression therapy to aid healing by effectively reducing swelling, edema, and pain:
  - Two reusable, anatomically shaped gel (water and glycerol gel) packs with open patella design for cooling
  - Integrated inflatable air chamber with pump for personalized compression

**APPLICATION INSTRUCTIONS - FIRST APPLICATION**

- Open all Velcro® straps.
- Apply and adjust the rods:
  - Insert the Velcro® pockets over the rods (fig.a1)
  - Attach the rods to the sides of the brace using Velcro® (fig.a2), ensuring the rods are aligned with the knee joint
  - If needed, shape the rods to match the patient’s leg anatomy
  - Adjust flexion-extension as prescribed: open the small cap, position the pins at the desired angle (fig.a3), then close the small cap
- Place the brace over the front of the knee, with the valve facing upward and the kneecap centered in the dedicated hole (fig.b).
- Fasten the brace using the Velcro® straps (fig.c). Use included extension straps for larger circumferences or added stability if necessary

**FOLLOW-UP APPLICATIONS**

- Open Velcro® straps
- Place the brace as described above (valve up, patella centered) (fig.b)
- Fasten using the Velcro® straps (fig.c)

**USE OF COLD GEL PACK**

- Remove the gel pack and place it in the freezer for at least two hours
- Attach the gel pack inside the brace using the dedicated Velcro® (fig.d)
- Apply the brace to the knee as previously described, always using a layer of fabric between the gel pack and the skin (e.g., the mesh sleeve included) to avoid direc contact wih cold surfaces. Do not use for more than 30 consecutive minutes.

**USE OF AIR CHAMBER**

- Screw the pump onto the appropriate valve on the brace(fig.e1).
- Turn the valve lever to “open” (fig.e2)
- Inflate the air chamber to the prescribed compression level (fig.e3) and close the valve by turning the lever to “close” (fig.e4).
- Unscrew the pump
- To partially deflate the air chamber, place the lever on “open” press the button on the pump (fig.e5). To fully deflate, remove the tube and leave the valve on “open” position.

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