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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

## Spinal orthosis for osteoporosis

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product. We recommend to not wear the device in the proximity of free frames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use. Contact an orthopaedic technician for the replacement of used parts.

### SELECTION/SIZE

Item	REF.26049	
Size	1	2
Distance T1-L5 cm	35/47	47/61
Colour	dark gray	

### MAINTENANCE

- Do not bleach  No chemical cleaning
  - Do not iron  Do not tumble-dry
  - Washing instructions: Remove the metal frame (see "Application", step 1). Hand wash in lukewarm water with neutral soap; rinse thoroughly. Allow to dry away from heat sources. Reinsert the metal frame (see "Application", step 1)
- Do not dispose of the device or any of its components in the environment.

### INDICATIONS

- Dorsolombalgia after healing of stable compression fractures of dorsolumbar vertebrae
- Stable compression fractures of dorsolumbar vertebrae caused by osteoporosis or metastatic tumors
- Hyperkyphosis of the thoracic spine, post-traumatic or caused by arthrosis
- Juvenile osteochondrosis
- Other indications upon medical prescription

### FEATURES AND MATERIALS

- Frame
  - Telescopic and moldable to the patient's anatomy
  - Made of radiolucent aluminum alloy and polymer
- Patented strap system
  - Abdominal grip with semi-corset
  - Unique tensioning system – pulling the terminals tensions all posterior and shoulder straps
  - Force-multiplying system that doubles traction at the shoulders – effective with minimal effort
- Dorsal-abdominal straps with low-friction belt; equipped with anatomical and ergonomic terminals
- Self-aligning return buckles
- Dorsal pressure point:
  - Adjustable with Velcro®
  - Elastic dorsal insert with bio-dynamic effect
- Shoulder straps with Velcro®-positionable padding
- Dorsal padding:
  - Made of bi-laminated PE with velvet and technical fabric; welded edges
  - Adjustable length
- Abdominal semi-corset for improved comfort and stability:
  - Anatomical design
  - Elastic part in bi-laminated MTP with technical fabric
  - Velcro® closure
  - Handles to facilitate closure
- Included:
  - Pair of additional ring pads
  - Pair of additional tubular pads

### APPLICATION

#### FIRST APPLICATION

- 1 Adjust the frame:
  - Release the two snap buttons (fig.A) at the lower base of the padding
  - Release the upper buckles (fig.B) by pushing them off the pin
  - Remove the frame from the padding (fig.C)
  - Place the frame on the patient's back to assess modifications (fig.D)
  - Adjust height by removing the pin (fig.E1), sliding the telescopic part to the desired length (fig.E2), and reinserting the pin (fig.E3)
  - Proceed with shaping using the edge of a table or, preferably, clamps (fig.F)
  - Reinsert the metal structure into the padding (fig.G), ensuring pins are at the top
  - Fasten the two buttons (fig.H)
  - Fasten the upper buckles (fig.I) until they click
- 2 Have the patient wear the brace by inserting arms into the shoulder straps (fig.J).
- 3 Abdominal adjustment:
  - Close the semi-corset with Velcro® at maximum length (fig.K)
  - Tighten abdominal straps (right and left), adjust length by trimming excess, and secure with Y-terminals (fig.L)
  - Ensure posterior structure is centered on the spine
- 4 Adjust shoulder straps using the Y-terminal (fig.M). Important: the attachment ring to the dorsolumbar strap must be centered on the midline (fig.N). Trim excess length.
- 5 Adjust padding so tubular part lies in the axillary cavity. Upper padding has small Velcro® for optimal fixation (fig.O).
- 6 Position dorsolumbar strap terminals beside abdominal pads (fig.P), check strap length under tension, and trim excess (fig.Q).
- 7 Verify overall brace position and strap tension, trimming further if necessary.

#### SUBSEQUENT APPLICATIONS

##### How to wear the brace:

- 1 Insert arms into shoulder straps (fig.R)
- 2 Close the semi-corset with Velcro®, ensuring central alignment (fig.S)
- 3 Release dorsolumbar strap terminals (pre-positioned at posterior end of abdominal straps), pull outward (fig.T), and fix with Velcro® at the anterior closure (fig.U)

##### To remove the brace:

- 1 Release dorsolumbar strap terminals (fig.V) and reposition them at posterior end of abdominal straps
- 2 Loosen shoulder straps by pulling forward (fig.W)
- 3 Release abdominal pads (fig.Y)
- 4 Remove the brace (fig.Z)

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