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
unicARPO 19+23

unicARPO 21+23

REF.26219 - Short-size wrist immobiliser, open, with rigid and malleable stays  
REF.26221 - Long-size wrist immobiliser, open, with rigid and malleable stays  
REF.26223 - Thumb accessory for UniCarpo braces

ORTHOSERVICE  
RO+TEN  
Take care feel better

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OSUK Ltd  
Units 47&50 North East BIC  
Wearfield  
Sunderland Enterprise Park East  
Sunderland SR5 2TA  
Tel +44 (0) 191 51 66 220  
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.  
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia  
Sede operativa e amministrativa:  
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia  
Tel. 0039 0362 18 51 200  
info@roplusten.com · www.roplusten.com

ORTHOSERVICE  
RO+TEN  
Take care feel better

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
## PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF.26219 - Short-size wrist immobiliser, open, with rigid and malleable stays

REF.26221 - Long-size wrist immobiliser, open, with rigid and malleable stays

REF.26223 - Thumb accessory for UniCarpo braces

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.

If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Please read the composition of the product on the inside label carefully; do NOT use the product if the patient has known intolerances or allergies to one or more components in the product.

We recommend to not wear the device in the proximity of free flames.

Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF.26219 • REF.26221 • REF.26223
Size	Universal
Length cm (UniCarpo 19)	22
Length cm (UniCarpo 21)	28
Colour	dark grey

Ambidextrous

### MAINTENANCE

- Do not bleach
- Do not iron
- Hand wash at max. 30° C
- The product is maintenance-free.
- Washing instructions: Remove the splint. Hand wash in lukewarm water with mild soap (at max. 30°C); rinse thoroughly. Dry away from heat sources.
- Do not dispose of the device or any of its components.

### INDICATIONS

#### UniCarpo19

- Wrist sprain
- Stable fracture of forearm radius and ulna
- Wrist tendinitis
- Arthritis and arthrosis
- After wrist fractures or plaster's early removal

#### UniCarpo21

- Treatment of stable fractures of meta-distal epiphysis of radius and ulna (deemed stable by orthopedist)
- Post-surgery after unstable distal fractures of radius and ulna
- Wrist sprain
- Tendinitis

#### UniCarpo21 + UniCarpo23

- Fractures of carpus (except scaphoid) whether treated surgically or conservatively
- Complications of rheumatoid arthritis
- Post-surgery after metacarpal fractures
- Radio-carpal arthrosis

#### UniCarpo19 + UniCarpo23

- Wrist sprain
- Wrist tendinitis
- Arthritis and arthrosis
- After wrist fractures or plaster's early removal
- Arthrosis of trapezio-metacarpal joint and metacarpophalangeal of thumb
- Damage to ulna section of thumb joint

#### UniCarpo21 + UniCarpo23

- A compound fracture of the first metacarpal
- A sprain of the radial ulnar collateral ligament of the metacarpophalangeal
- A suspected scaphoid fracture

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

#### UniCarpo19, 21

- Base made of doublelaminated MTP with soft velour fabric
- Universal and ambidextrous design: one product for all needs
- Rigid palmar aluminum stay, pre-shaped and moldable according to medical instructions
- Dorsal reinforcement stays in aluminum
- Closure with Velcro® straps and return rings; length adjustable as needed thanks to Yshaped strap ends
- Full dorsal opening

#### UniCarpo23

- Thumb support attachable with Velcro® to UniCarpo19 and UniCarpo21 braces
- Base made of doublelaminated MTP with soft velour fabric
- Aluminum thumb stay, reshaped and moldable according to medical instructions
- Velcro® closure

### APPLICATION OF UNICARPO19 AND UNICARPO21

#### 1 ONLY FOR THE FIRST APPLICATION

- Shape the rigid palmar stay (Fig. A) to fit the patient's wrist after removing it from the orthosis. Reinsert the stay into the dedicated pocket.
  - Shape the dorsal stays as well.
- Loosen all the Velcro® straps (temporarily fastening them, onto themselves) and if necessary fully open the brace.
- Position the brace so that the rigid palmar stay is on the palm side of the hand and the thumb fits into the dedicated opening. (Fig. B)
- Secure the brace:
  - Fasten the central strap with Velcro®. (Fig. C1)
  - Fasten the remaining straps and secure with Velcro®: the hand strap (Fig. C2) and the forearm one. (Fig. C3)
  - (ONLY ref.26221) Tighten and fasten the proximal forearm strap with Velcro® as well. (Fig. D)
- If the straps are too long, they can be shortened by cutting them to the desired length and reattaching the corresponding Yshaped end.

#### APPLICATION UNICARPO23

##### 1 ONLY FOR THE FIRST APPLICATION

- Attach the UniCarpo23 thumb support to the UniCarpo19 or UniCarpo21 brace using the Velcro®. Center it over the opening (right or left, depending on the affected side) of the base wrist cuff. (Fig. E)
  - Adjust the thumb opening using the dedicated Velcro® strap. (Fig. F)
  - If necessary, shape the thumb stay.
- Apply the brace as described in steps 1-5.

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