

cervistrong

Reinforced neck brace

REF. 20304

Indications

- Mild or moderate distortions of the cervical spine (whiplash)
- Cervical vertigo
- Cervical osteoarthritis
- Moderate discopathies with associated cervicobrachial pain syndrome

RADIO
LUCENTANATOMICAL
FIT

Code	REF. 20304			
Size	S	M	L	XL
Neck circumference cm	30/35	36/41	42/47	48/53
Height cm	6-8-10-12	6-8-10-12	6-8-10-12	6-8-10-12
Colour	grey			



orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

Siedziba w Polsce: Orthoservice Polska Sp. z o. o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10
info@orthoservice.pl · www.orthoservice.pl

Sede italiana: RO+TEN S.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland

cervistrong

Reinforced neck brace



**ORTHO SERVICE
RO+TEN**

Take care feel better



**ORTHO SERVICE
RO+TEN**

Take care feel better

A painful neck

The cervical region is the least protected part of the spine, so it may suffer from **distortions** (strained muscles and tendons) or **osteoarticular diseases** (arthritis and discopathies) which can be particularly debilitating.

An innovative solution

Even the most ordinary movements, such as turning the head or slightly bending the neck, trigger pain that can also radiate to the shoulders and arms, causing an intense feeling or stiffness and sometimes even dizziness.



Limiting any movement as much as possible is the treatment of choice, but until now that choice was between a rigid neck brace, effective but uncomfortable, and a soft neck brace, comfortable but providing less support. The innovatively designed **Cervistrong combines levels of efficacy and comfort that have never been seen before.**

THE SUPPORT OF A SEMIRIGID NECK BRACE

A **special stiff polyethylene core** with polymer-reinforced edges immobilises the cervical spine helping it to support the neck as it normally would. It feels imperceptible against the skin.


THE COMFORT OF A FOAM NECK BRACE

The soft foam rubber design is already shaped and **adapts to fit the wearer's neck** by following the curves of the jaw and the nape of the neck.

RELIABLE, COMFORTABLE AND PRACTICAL

The brace is covered with a **soft jersey fabric** that feels soft and comfortable against the skin. It fastens at the back with Velcro®, making the brace quick and easy to put on.

Speeds up healing and relieves pain



Sturdy polyethylene core with polymer edges.

The soft foam rubber design is anatomically shaped.

Jersey covering and rear Velcro® fastening.