

PRODUCT LINE **sat**

The evolution of success
is an excellence



**ORTHOSERVICE
RO+TEN**

Take care feel better

The great new additions of the most loved and imitated product line



THE DESIGN OF THE ELASTIC STRAPS

The **elastic fabric** has a **double parallel band** that ensures traction in the center and softness on the edges. Hypoallergenic (latex-free), it does not unravel or sag. The traditional traction system allows for **double the force during closure** compared to the same amount of force exerted.



A FASHIONABLE TOUCH OF TEXTURE

The modern and dynamic texture of the fabric finally gives the **aesthetic pleasure** that promotes self-confidence.



CONTAINMENT, SEALING AND FIT OF THE FABRIC

Thanks to the new **reinforced fabric, 98% cotton**, the orthopedic effectiveness, comfort, and breathability, have reached a unique and inimitable quality.



THE COMFORT OF THE LOWER SHAPING

The profile of the brace at the level of the leg joint has been **shaped** to perfectly follow the lines of the body.



2 PADS IN EACH PACK

The abdominal and back pads for stabilization are ready for use.



A product line ready for the challenges of the future



THE EASY GRIP OF THE ELASTIC TIE RODS

The **soft polymeric material** of the tie-rods facilitates the grip both while wearing and removing the brace.



THE LOW PROFILE AND PRACTICALITY OF THE FRONT CLOSURE

The closure has been redesigned so that the edges do not overlap. Traction is aided by the new **thumb grip loop**.



THE HYPOALLERGENIC PADDING OF THE SHOULDER STRAPS

Redesigned and made of **latex-free fabric**, the padding ensures optimal comfort.
























REGULAR FIT AND TAILORED FIT, THE NEW INTUITIVE SIZES

The new size definition makes the choice even easier: **Regular Fit** indicates a brace with a standard fit; **Tailored Fit** is suitable for those with a narrower waist.



The entire product line at a glance

		TAILORED FIT						REGULAR FIT							
		Size	XS	S	M	L	XL	XXL	Size	XS	S	M	L	XL	XXL
		Pelvic circ. cm	75/85	85/95	95/105	105/115	115/125	120/135	Pelvic circ. cm	70/80	80/90	90/100	100/110	110/120	120/130
 <p>Indications</p> <ul style="list-style-type: none"> • Lumbago • Sciatic lumbalgia and lumbar cruralgia • Spondylarthrosis and lumbar discopathies • Slight traumas to the lumbosacral column and paravertebral lumbar contractures • Post-operative recovery • After-effects of stable compression fractures of lumbar vertebrae • Treatment of fracture of transverse processes 	<p>sat 1323 low</p>  <p>Back height 33 cm</p>  <p>Front height 21 cm</p>							<p>sat 1355 low</p>  <p>Back height 32 cm</p>  <p>Front height 20 cm</p>							
	 <p>Indications</p> <ul style="list-style-type: none"> • Lumbalgia • Sciatic lumbalgia and lumbar cruralgia • Spondylarthrosis and lumbar discopathies • Slight traumas to the lumbosacral column and paravertebral lumbar contractures • Post-operation period • Results of vertebral lumbar relaxations • Treatment of fracture of transverse processes • Degenerative lumbar scoliosis • Spondylosis with or without vertebral listhe 	<p>sat 1322 high</p>  <p>Back height 38 cm</p>  <p>Front height 23 cm</p>	<p>sat 1322 SH high front lowered</p>  <p>Back height 38 cm</p>  <p>Front height 19 cm</p>							<p>sat 1354 high</p>  <p>Back height 35 cm</p>  <p>Front height 23 cm</p>					
		<p>sat 1321 prominent abdomen</p>  <p>Back height 38 cm</p>  <p>Front height 31 cm</p>							<p>sat 1353 prominent abdomen</p>  <p>Back height 35 cm</p>  <p>Front height 31 cm</p>						
	 <p>Indications</p> <ul style="list-style-type: none"> • Back pain due to arthrosis or kyphosis • After-effects of fractures of dorsal vertebrae and/ or upper lumbar vertebrae • After-effects of stable compression fractures of the dorsal spine due to osteoporosis or metastasis 	<p>sat 1336 padded</p>  <p>Back height 55 cm</p>  <p>Front height 21 cm</p>							<p>sat 1335 padded</p>  <p>Back height 55 cm</p>  <p>Front height 23 cm</p>						

Low orthosis

with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

REF. 1323 tailored fit

REF. 1355 regular fit

Option

9ST012/L optional rigid stays cm 26



Tailored fit
Front height
21 cm
Back height
33 cm

Regular fit
Front height
20 cm
Back height
32 cm

High and high front lowered orthosis

with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

REF. 1322 tailored fit

REF. 1322SH tailored fit
front lowered

REF. 1354 regular fit

Option

9ST012/L optional
rigid stays cm



Tailored fit
Front height
23 cm
Back height
38 cm

**Tailored fit,
front lowered**
Front height
19 cm
Back height
38 cm

Regular fit
Front height
23 cm
Back height
35 cm



Orthosis for prominent abdomen

for stabilising and supporting the lumbar rachis, with rigid and thermoformable pads

REF. 1321 tailored fit

REF. 1353 regular fit

Option

9ST012/L optional rigid stays cm 30

Tailored fit
Front height 31 cm
Back height 38 cm



LATERAL OPENING

To make the brace for prominent abdomen issues easy and simple to put on, a side opening has been provided.



Regular fit

Front height 31 cm
Back height 35 cm



Padded orthosis

for stabilising and supporting the dorso-lumbar rachis, with rigid and thermoformable pads

REF. 1336 tailored fit

REF. 1335 regular fit

Option

9ST012/L optional rigid stays cm 38



Tailored fit
Front height 21 cm
Back height 55 cm



Regular fit

Front height 23 cm
Back height 35 cm





ORTHOSERVICE
RO+TEN

Take care feel better

LINEA sat

Orthoses with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

Available in different versions:

low • high • prominent abdomen • padded

9DCDP001581TA01.0322

TAILORED FIT

Size	XS	S	M	L	XL	XXL
Pelvic circ. cm	75/85	85/95	95/105	105/115	115/125	120/135

Item REF. 1323

Front height cm	21
Back height cm	33

Item REF. 1322

Front height cm	23
Back height cm	38

Item REF. 1322 SH

Front height cm	19
Back height cm	38

Item REF. 1321

Front height cm	31
Back height cm	38

Item REF. 1336

Front height cm	21
Back height cm	55

REGULAR FIT

Size	XS	S	M	L	XL	XXL
Pelvic circ. cm	70/80	80/90	90/100	100/110	110/120	120/130

Item REF. 1355

Front height cm	20
Back height cm	32

Item REF. 1354

Front height cm	23
Back height cm	35

Item REF. 1353

Front height cm	31
Back height cm	35

Item REF. 1335

Front height cm	23
Back height cm	55



ANATOMICAL
FIT



BREATHABLE



CUSTOMIZABLE
FIT



EASY
PULLING



SINGLE PATIENT
MULTIPLE USE



RIGID
THERMOFORMABLE
PAD

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

Siedziba w Polsce: Orthoservice Polska Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10
info@orthoservice.pl · www.orthoservice.pl

Sede italiana: RO+TEN S.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:

Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



orthoservice.com