

shouldfix light

Orthosis for shoulder immobilization

| | |
|--------|----------------|
| Code | REF.20413 |
| Size | universal size |
| Colour | black |

fits right and left

Indications

- Postoperative immobilisation following subacromial decompression, surgical removal of calcifications, sutures for the rotator cuff, revision surgery, arthrolysis, osteosynthesis, reverse shoulder replacement, and many other conditions
- All clinical cases that require early functional exercise (free movement of the elbow and wrist)
- Contusions
- Distortions

Contraindications

- Currently no known



shouldfixlight

Orthosis for shoulder immobilization



www.orthoservice.com

UNI EN ISO 13485:2016
UNI EN ISO 9001:2015



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

Sede italiana: RO+TEN s.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



**ORTHOSERVICE
RO+TEN**

Take care feel better



**ORTHOSERVICE
RO+TEN**

Take care feel better



shouldfix light

Reduce pain and promote healing

Shouldfix light is used in the conservative and postoperative treatment of a range of shoulder conditions.

This brace has been designed to immobilise the shoulder in the most comfortable way possible yet without unduly restricting movement.

According to the latest theory about functional exercise, it has been shown that conventional



immobilisation with internal rotation of the arm can cause pain, tension and stiffness: this is why the basic configuration does not provide for fixing the forearm and wrist

Shouldfix light is the optimum solution in all clinical cases that require early functional exercise, such as distortion or contusion of the shoulder.



ELBOW AND WRIST MOBILITY

The brace comprises an abdominal strap and an armband with Velcro®.

Made with nonslip inserts, the strap can be adapted to fit the patient and is easy to adjust. The armband is positioned above the elbow, and immobilises the arm while still allowing the elbow and wrist to move freely.

An immobilised shoulder heals better with freedom of movement



OPTIONAL FOREARM RESTRICTION

The patient's forearm and wrist can also be secured to the abdominal strap with easily positionable Velcro® bands.



The entire brace is made with a hassle-free fine, lightweight fabric that can even be worn against the skin without any problems, guaranteeing a comfortable fit and practical design. The brace can be easily washed.