

shouldfix 3

Orthosis for shoulder immobilization

Code	REF.20414			
Size	S	M	L	XL
Forearm length* cm	28/34	34/42	42/48	48/58
Colour	black			

Fits right and left
 *measurement taken from the elbow to the end of the closed fist

Indications

- Contusions
- Distortions
- Anterior dislocations
- Simple fractures
- Postoperative immobilisation following shoulder arthroscopy, endoprosthesis or treatment of a fracture

Contraindications

- Currently no known



shouldfix 3

Orthosis for shoulder immobilization



www.orthoservice.com

UNI EN ISO 13485:2016
 UNI EN ISO 9001:2015



Headquarter: ORTHOSERVICE AG
 Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
 Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
 info@orthoservice.com · www.orthoservice.com
Niederlassung Deutschland: Orthoservice Deutschland GmbH
 Flugstraße 8 · D-76532 Baden-Baden · Deutschland
 Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
 info@orthoservice.de · www.orthoservice.de
Sede italiana: RO+TEN s.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
 Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
 Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
 info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):
 Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



**ORTHOSERVICE
 RO+TEN**

Take care feel better



**ORTHOSERVICE
 RO+TEN**

Take care feel better



shouldfix3 is the ideal solution for immobilising and stabilising a shoulder joint **in the conventional internal rotation position**

When healing depends on **precision and comfort**

The brace must always be correctly applied whenever a shoulder has been immobilised in internal rotation. It is therefore of the utmost importance that the patient is capable of donning and doffing the brace with ease.



THE CORRECT FITTING GUARANTEED

To make using the brace as easy as possible for patients and medical staff, we have colour-coded the Velcro® fastenings, rings and straps: black for the shoulder strap and white for the strap around the waist. It is a foolproof solution.