

shouldercross

Functional shoulder
brace

NEW



ORTHOSERVICE
RO+TEN

Take care feel better

An ally for the unstable shoulder: **Shouldercross**

SHOULDER INSTABILITY CAN BE:

- **Idiopathic:** congenital excess laxity and elasticity of the ligamentous structures and capsule
- **Microtraumatic:** common in athletes at the extremes caused by repetitive athletic movements
- **Post-traumatic:** complete or partial loss of the joint relationship between the humeral head and the glenoid.

The new **Shouldercross**, version is the innovative therapeutic alternative

For preventive or post-traumatic phases

Wearing **Shouldercross** for 2 to 6 hours during “at-risk” sports and work activities helps stabilize the joint: it limits dislocation and subluxation episodes, restoring confidence in the patient.

For painful loose shoulders

Using **Shouldercross** for 2 to 4 months, combined with proper rehabilitation, reduces recovery time and the percentage of necessary stabilization surgeries.

In post-operative phases

Shouldercross allows for a much more tolerable recovery and rehabilitation phase: it assists without immobilizing.

In all cases

The use of **Shouldercross** is strongly recommended while the patient perceives benefit.

NEW ADJUSTABLE ARM GRIP

The arm grip can be adjusted with a **practical Velcro® closure** to perfectly fit the limb.

NEW



CROSS-REINFORCED ELASTIC STRAPS

Limit the desired movement
according to medical prescription:

Position 1 or

Position 2



**CONTAINMENT ALSO ON
THE POSTERIOR PART**

BREATHSKILL FABRIC STRUCTURE

Soft and breathable thanks
to its **special 3D design**.

The external mesh weave
provides strength.



Innovative and
functional, it stabilizes
the shoulder without
immobilizing it.



ADJUSTABLE AND PADDED CONTRALATERAL FIXING STRAP

Easy to close with a **spring balancer**.
It stabilizes the orthosis and
prevents migration.
The strap can be cut to size,
ensuring the padding is placed
under the armpit.



shouldercross

Functional shoulder brace

NEW

REF.20726 R

REF.20726 L

Size	S	M	L	XL	XXL
Circ. biceps cm	22/26	26/30	30/35	35/40	40/45
Colour	grey				

Specify right or left



RADIO
LUCENT



LOW
PROFILE



ANATOMICAL
FIT



BREATHABLE



3D FABRIC



SINGLE PATIENT
MULTIPLE USE



Indications

- Conservative treatment of instability after both front and back glenohumeral dislocations
- Protection during high-risk sports or work activities, while awaiting surgery or after stabilization surgery
- Use for shoulder rehabilitation treatment of painful hyperlaxity or functional instability
- Useful after repair of large rotator cuff tears or for conservative management of massive rotator cuff tears
- Useful in the rehabilitation management of scapular dyskinesia.

Contraindications

- Currently none.



www.orthoservice.com



UNICI EN ISO 15488:2021
UNI EN ISO 9001:2015



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**

Take care feel better