

# elbo 2.0

All elbow joint conditions need to be treated with the utmost care irrespective of whether they have been caused by injury. **An immobilising brace** is the optimal solution during immediate **post-operative** care to temporarily restrict **flexion, extension, pronation and supination**, and also during rehabilitation when joint stability plays an essential role in aiding tissue repair and enabling the joint to function correctly again. In both cases the brace must be **effective and comfortable** as well as adjustable and easily managed by the patient. This the only way to ensure a protected, **easy and safe path to recovery.**



9003P00116EN003.1220

## elbo 2.0

Telescopic elbow brace with R.O.M.

Code	REF.B1017
Size	one size
Color	black

specify right and left

### Indications

- Compound fracture of the elbow, distal humerus, and proximal radius and ulna.
- Postoperative immobilisation following surgery on tendons and ligaments.
- Beneficial during rehabilitation following surgery and conservative treatment of elbow fractures.
- Subluxation or dislocation of the elbow.
- Chronic elbow injuries.
- Hyperextension of the elbow.
- Post-trauma treatment.

### Contraindications

- There are no known contraindications at this time.



Option: M.B1018  
NEW OPTION  
NEUTRAL HAND SUPPORT



# elbo 2.0

**Telescopic  
elbow brace  
with R.O.M.**



NEW



www.orthoservice.com



**Headquarter:** ORTHOSERVICE AG  
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland  
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89  
info@orthoservice.com · www.orthoservice.com  
**Niederlassung Deutschland:** Orthoservice Deutschland GmbH  
Flugstraße 8 · D-76532 Baden-Baden · Deutschland  
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13  
info@orthoservice.de · www.orthoservice.de  
**Sede italiana:** RO+TEN s.r.l.  
**Sede legale:** Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia  
**Sede operativa e amministrativa:**  
Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia  
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34  
info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):  
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



**ORTHOSERVICE  
RO+TEN**

Take care feel better



**ORTHOSERVICE  
RO+TEN**

Take care feel better



# NO SCISSORS REQUIRED! ELBOW BRACES HAVE NEVER BEEN THIS EASILY ADAPTABLE!

NEW

## A NEW TELESCOPIC STAY FOR PERFECTLY CUSTOMISED ADJUSTMENTS

The lightweight aluminium alloy lateral stay is telescopic and adjustable in two points: above and below the elbow. Just gently pull the lever on the cuffs and slide them along the stay to the desired position to suit the size of the joint. A click indicates that it has locked into the desired hole. The system is safe and cannot be accidentally unlocked. Once the brace is carefully fitted, no more worries about cutting straps to fit.

## A SIMPLIFIED NEW HINGE WITH AN ANTI- RELEASE SAFETY LOCK

Made from polycarbonate, the hinge for adjusting the flexion and extension is ever-so easy to work: just pull the sliders outwards (white for flexion and orange for extension) and turn the disc until you reach the desired degree of flexion/extension and then push them back in. The system is also very safe; once the setting has been locked in, it cannot be accidentally changed by the patient if the brace knocks against something. Adjustable extension from 0° to 90°. Adjustable flexion from 0° to 120°.

## NEW OPTION NEUTRAL HAND SUPPORT

## THERMOPLASTIC CUFFS

Two fixed and two mobile. The cuffs have been specially designed to accommodate the patient's arm and forearm and immobilise them without applying any pressure. The length of the straps can be easily adjusted: just remove the Y-tab, cut off any excess, and reattach the tabs to the ends.



## THERMOFORMED PADDING

Made from bonded MTP, the thermoformed padding guarantees a comfortable fit and prevents any annoying irritations from occurring anywhere on the skin.

## ADJUSTABLE SHOULDER STRAP FOR SUPPORT

The supportive shoulder strap is incredibly easy to adjust and wear: just pass the Velcro end through the ring provided on the wrist, fold it back on itself and attach it to the Velcro®. Once it has been adjusted to the right length, the patient can don the shoulder strap without any assistance just by opening and closing the clip at the opposite end which attaches to the same ring. Made from soft velour, it features padded shoulder protection.



elbo 2.0