

# polFIT LINE



## polFIT 19

REF.P1119 R/L

Short-size wrist immobilizer, open, with a rigid and malleable stay

Code	REF.P1119 R/L · REF.P1117 R/L				
Size	XS	S	M	L	XL
Wrist circumf. cm	11/13	13/16	16/18	18/20	20/23
Length cm	20	20	20	20	20
Colour	black				

specify right or left



## polFIT 17

REF.P1117 R/L

Short-size wrist immobilizer with thumb, open, with rigid and malleable stays



## polFIT 21

REF.P1121 R/L

Long-size wrist immobilizer, open, with a rigid and malleable stay

Code	REF.P1121 R/L · REF.P1122 R/L			
Size	S	M	L	XL
Wrist circumf. cm	13/16	16/18	18/20	20/23
Length cm	26	26	26	26
Colour	black			

specify right or left



## polFIT 22

REF.P1122 R/L

Long-size wrist immobilizer with thumb support, open, with rigid and malleable stays



www.orthoservice.com



**Headquarter: ORTHOSERVICE AG**  
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland  
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89  
info@orthoservice.com · www.orthoservice.com



**ORTHOSERVICE DEUTSCHLAND GmbH**  
Flugstraße 8 · D-76532 Baden-Baden · Deutschland  
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13  
info@orthoservice.de · www.orthoservice.de

**ORTHOSERVICE POLSKA Sp. z o.o.**  
ul. Warszawska 416a · 42-209 Częstochowa · Polska  
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785  
info@orthoservice.pl · www.orthoservice.pl

**OSUK Ltd**  
Units 47&50 North East BIC  
Wearfield  
Sunderland Enterprise Park East  
Sunderland SR5 2TA  
Tel +44 (0) 191 51 66 220  
info@orthoservice.uk · www.orthoservice.uk

**RO+TEN S.r.l.**  
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia  
Sede operativa e amministrativa:  
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia  
Tel. 0039 0362 18 51 200 · Fax 0039 0362 18 58 065  
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE  
RO+TEN**

Take care feel better

# polFIT LINE

## A complete line of wrist immobilizers



**ORTHOSERVICE  
RO+TEN**

Take care feel better

**SHORT-SIZE**  
REF.P1119 R/L



**SHORT-SIZE WITH THUMB**  
REF.P1117 R/L



**LONG-SIZE**  
REF.P1121 R/L



**LONG-SIZE WITH THUMB**  
REF.P1122 R/L



LENGTH 20 CM

LENGTH 26 CM

## The characteristics of the **polFIT** line



1 Full dorsal opening



2 Velcro® fastening system



3 3D breathable fabric for structure



4 Aluminum palm splint, preformed and removable on medical prescription



5 Reinforcing dorsal stays



6 **ONLY PoIFIT 17 and PoIFIT 22** Aluminum thumb splint, preformed and mouldable on medical prescription



BREATHABLE



EASY FASTENING



3D FABRIC