

POIFIT LINE

**A complete line
of wrist
immobilizers**



**ORTHOSERVICE
RO+TEN**

Take care feel better

SHORT-SIZE
REF.P1119 R/L



SHORT-SIZE WITH THUMB
REF.P1117 R/L



LONG-SIZE
REF.P1121 R/L



LONG-SIZE WITH THUMB
REF.P1122 R/L



LENGTH 20 CM

LENGTH 26 CM

The characteristics of the **polFIT** line



1 Full dorsal opening



2 Velcro® fastening system



3 3D breathable fabric for structure



4 Aluminum palm splint, preformed and removable on medical prescription



5 Reinforcing dorsal stays



6 **ONLY PoIFIT 17 and PoIFIT 22** Aluminum thumb splint, preformed and mouldable on medical prescription



BREATHABLE



EASY FASTENING



3D FABRIC

polFIT LINE



polFIT 19

REF.P1119 R/L

Short-size wrist immobilizer, open, with a rigid and malleable stay



polFIT 17

REF.P1117 R/L

Short-size wrist immobilizer with thumb, open, with rigid and malleable stays

Code

Size

Wrist circumf. cm

Length cm

Colour

REF.P1119 R/L · REF.P1117 R/L

S

M

L

XL

13/16

16/18

18/20

20/23

20

20

20

20

black

specify right or left



polFIT 21

REF.P1121 R/L

Long-size wrist immobilizer, open, with a rigid and malleable stay



polFIT 22

REF.P1122 R/L

Long-size wrist immobilizer with thumb support, open, with rigid and malleable stays

Code

Size

Wrist circumf. cm

Length cm

Colour

REF.P1121 R/L · REF.P1122 R/L

S

M

L

XL

13/16

16/18

18/20

20/23

26

26

26

26

black

specify right or left



SINGLE PATIENT
MULTIPLE USE



BREATHABLE



3D FABRIC



EASY
FASTENING



RADIO
LUCENT



www.orthoservice.com

UNI EN ISO 13485:2016
UNI EN ISO 9001:2015



Organismo accreditato da ACCREDIA
Body accredited by ACCREDIA



ORTHOSERVICE
RO+TEN

Eco
attentive

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

Siedziba w Polsce: ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa Polska
Tel.: +48 (0) 34 340 13 10 · (NIP): 9492246785
www.orthoservice.pl · info@orthoservice.pl

Sede italiana: RO+TEN S.r.l.

Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia

Sede operativa e amministrativa:

Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento
(art. 2497bis CC): Orthoservice AG (CH)
6830 Chiasso (TI) · Switzerland



**ORTHOSERVICE
RO+TEN**

Take care feel better