


polFIT LINE



polFIT 17

REF.P1117 R/L

Short-size wrist immobilizer with thumb support, open, with rigid and malleable stays

| | | | | | |
|-------------------|-------------------------------|-------|-------|-------|-------|
| Code | REF.P1117 R/L · REF.P1119 R/L | | | | |
| Size | XS | S | M | L | XL |
| Wrist circumf. cm | 11/13 | 13/16 | 16/18 | 18/20 | 20/23 |
| Length cm | 20 | 20 | 20 | 20 | 20 |
| Colour | black | | | | |

specify right or left



polFIT 19

REF.P1119 R/L

Short-size wrist immobilizer, open, with rigid and malleable stays



polFIT 21

REF.P1121 R/L

Long-size wrist immobilizer, open, with rigid and malleable stays

| | | | | |
|-------------------|-------------------------------|-------|-------|-------|
| Code | REF.P1121 R/L · REF.P1122 R/L | | | |
| Size | S | M | L | XL |
| Wrist circumf. cm | 13/16 | 16/18 | 18/20 | 20/23 |
| Length cm | 26 | 26 | 26 | 26 |
| Colour | black | | | |

specify right or left



polFIT 22

REF.P1122 R/L

Long-size wrist immobilizer with thumb support, open, with rigid and malleable stays



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**
Take care feel better

polFIT LINE


A complete line of wrist immobilizers



**ORTHOSERVICE
RO+TEN**
Take care feel better

SHORT-SIZE WITH THUMB

REF.P1117 R/L - Short-size wrist immobilizer with thumb support, open, with rigid and malleable stays




polFIT 17

LENGTH 20 CM

SHORT-SIZE

REF.P1119 R/L - Short size wrist immobilizer, open; with rigid and malleable stays



polFIT 19

LONG-SIZE

REF.P1121 R/L - Long size wrist immobilizer, open; with rigid and malleable stays




polFIT 21

LENGTH 26 CM

LONG-SIZE WITH THUMB

REF.P1122 R/L - Long size wrist immobilizer with thumb support, open; with rigid and malleable stays



polFIT 22

The characteristics of the polFIT line



1 Full dorsal opening



2 Velcro® fastening system



3 3D breathable fabric for structure



4 Aluminum palm splint, preformed and removable on medical prescription



5 Reinforcing dorsal stays in aluminum, mouldable



6 **ONLY PolFIT 17 and PolFIT 22** Aluminum thumb splint, preformed and mouldable on medical prescription



BREATHABLE



EASY FASTENING



3D FABRIC