

ONE BRACE, MULTIPLE APPLICATIONS:

post-arthroscopy care	▶ oedema prevention, cryotherapy and movement control
conservative treatment of hip dislocations	▶ movement control (flexion and extension and especially postoperative abduction/adduction)
post-surgical care (also for replacements or revision surgery)	
arthrosis of the hip joint	▶ pain relief
coxalgia or inoperable joint misalignment	

9DCDP00027EN003_0421



REF. 91007
Adjustable rod for Hipocross

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de
ORTHOERVICE POLSKA sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa Polska
Tel.: +48 (0) 34 340 13 10
www.orthoservice.pl · info@orthoservice.pl
Sede italiana: RO+TEN s.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com
Società soggetta a Direzione e Coordinamento (art. 2497bis CC);
Orthoservice AG (CH) 6830 Chiasso (TI) · Switzerland

Hipocross

Brace for hip and coxalgia

Item	REF. 91006 (Hipocross)			
Size	S	M	L	XL
Pelvis circ. cm	72/93	93/116	116/140	140/165
Waist circ. cm	65/84	84/109	109/140	140/170
Proximal thigh circ. cm	40/54	50/65	60/76	70/87
Colour	black			
Specify right or left				

Item	REF. 91007 (Adjustable rod for Hipocross)
Size	one size
Colour	black
Specify right or left	

INDICATIONS

- Post-arthroscopy care.
- Post-surgical care (compatible with the degree of stabilization required).
- Conservative treatment of hip dislocations.
- Arthrosis of the hip joint
- Coxalgia or inoperable joint misalignment

CONTRAINDICATIONS

Currently no known



ORTHOERVICE
RO+TEN
Take care feel better

HIPOCROSS, THE X THAT MAKES THE DIFFERENCE

Light weight

Stability

Comfort

Safety



ORTHOERVICE
RO+TEN
Take care feel better



RIGID BAR

A rigid bar that locks the joint in flexion and neutral abduction during the initial days following a dislocation.

ADJUSTABLE ROD, THE ADDED BONUS FOR ADDED SAFETY

The adjustable rod enables you to regulate:

- adduction by adjusting the hinge (positioned level with the trochanter) with a screwdriver and assisting the leg into the desired degree of adduction
- flexion and extension by setting the graduated hinge to the desired degree between 0° and 120°.



ATTENTION TO DETAIL

CRYOTHERAPY GEL AND COMPRESSION FABRIC

Thanks to compression fabric and gel that can be chilled in the freezer and placed in the special pocket over the joint, **Hipocross** aids the reabsorption of oedema.



STRAPS DESIGNED WITH FOUR OVERLAPPING STRETCHY BANDS

We have specially designed two sturdy straps with four overlapping stretchy bands that provide four times more oppositional force than an individual elasticated band so that they **cross** at the side and provide dynamic adduction control.

