



malleoLITE

ANKLE SUPPORT WITH ELASTICATED FIGURE-8 STRAPS

REF. 33080



AS LIGHT AS
A FEATHER, AS
SUPPORTIVE AS
A FIGURE-8
BANDAGE



ORTHOERVICE ORTHOPEDIC PRODUCTS



Made with lightweight Sensitive® fabric which is breathable, ladder-proof and antistatic.

The fully open design is ever-so easy to put on without having to flex the foot.

Anatomically designed with soft, fine finishes.

Elasticated mille-rayé figure-8 straps provide a valid alternative to a bandage.

The fine, lightweight design slips easily into a shoe.

The rear fastening Velcro® adjusts to fit around the ankle.

malleoLITE

ANKLE SUPPORT WITH ELASTICATED FIGURE-8 STRAPS

REF. 33080

REF. 33080

Size

S

M

L

Circ. of ankle from base of heel across the malleoli cm

24-28

28-33

33-38

Colour: black / Fits right and left



Headquarter: ORTHOSERVICE AG

Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 91 8220088 · Fax 0041 91 8220089

info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH

Flugstraße 8 · D-76532 Baden-Baden · Deutschland

Tel. 0049 (0)7221 9913911 · Fax 0049 (0)7221 9913913

info@orthoservice.de · www.orthoservice.de

Sede italiana: RO+TEN s.r.l.

Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia

Sede operativa e amministrativa:

Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia

Tel. 0039 039 6014094 · Fax 0039 039 6014234

info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento [art. 2497bis CC]:

Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland

INDICATIONS

- Grade 1 ankle sprains.
- Tendinopathy of the instep and peroneal tendons.
- Mild joint pain.



ORTHOSERVICE
R+TEN