

ORTHO 14

**Supporting
the next generation**



**ORTHOSERVICE
RO+TEN**

Take care feel better



PEDIATRIC LINE **ORTHO 14**

ORTHO 14-100



Soft foam cervical collar

- **CORE MADE** of semi-rigid foam
- **TUBULAR** jersey cover
- **VELCRO® FASTENING** at back of neck
- **X-RAY TRANSPARENT**

INDICATIONS

- Slight cervical spine distortions (whiplash)
- Stiff necks
- Neck-induced vertigo

Item	REF.OFT100
Size	universal
Height cm	6,5
Circonf. cm	20/30
Colour	beige



Bivalve cervical collar with tracheal hole

- Non-toxic PLASTAZOTE® **STRUCTURE**
- **CONSISTS OF TWO ENVELOPING VALVES** with Velcro® closure
- **SHAPED** for perfect fit at chin, neck and shoulders
- Central front and back rigid **SUPPORT PLATE**
- **HOLES** for ventilation
- **RADIOLUCENT**
- Trachea **HOLE**
- **PREDISPOSITION** for attachment of the thoracic stabilizer

INDICATIONS

- Medium severity to severe distortion of the cervical spine (whiplash)
- Multiple discopathies of the cervical spine
- Cervical spine traumas
- Pathological instability as a result of cervical spine metastases
- Support of post-surgical spinal surgery care

Item	REF.OFT105I	REF.OFT105P
Size	Kids 0 – 3 years old	Kids 4 – 9 years old
Colour	beige	



Immobiliser for arm and shoulder, with elbow support

- **SHOULDER STRAPS** with Velcro® fastenings
- Soft padded fabric **IMMOBILISER SWATHE**
- Elbow **SUPPORT**
- **SOFT PADDED** fabric immobiliser swathe

INDICATIONS

- Compound humerus fracture
- Glenohumeral subluxation
- Shoulder sprain
- Post-surgical care of humeral fractures and interventions on humeral subluxations
- Postoperative care of the long head of the biceps tendon

Item	REF.OFT305
Size	universal
Length bag cm	23
Height cm	22
Colour	black

Fits right and left



Arm sling

- Velcro® **FASTENING**
- **SHOULDER STRAP** with neck protection cell
- **STRUCTURE** made of cotton blended fabric

INDICATIONS

- Forearm and wrist support aid
- As support in cases of cast limb
- Useful if resting of the acromioclavicular (after a subluxation grade 1) or glenohumeral joint is required
- Useful in case of immobilisation of the upper arm with the elbow in flexed position

Item	REF.OFT310
Size	universal
Length bag cm	24
Colour	blue

Fits right and left



Clavicle splint

- Made of gauze fabric **WITH PADDING**
- **APPLIED** similarly to the figure of eight bandage with Velcro® closure at the back

INDICATIONS

- Clavicle fractures
- Subluxation of the acromionclavicular joint

Item	REF.OFT315
Size	universal
Colour	white

Fits right and left



Wrist immobiliser

- Aluminium palm **SPLINT**, preformed and pliable (removable)
- Reinforcing dorsal **STAYS**
- **FASTENING** with Velcro® straps and buckles
- Double-lined foam **DESIGN**

INDICATIONS

- Wrist distortion
- Compound fracture of the forearm of the radius and ulna extremities
- Wrist tendonitis
- Consequences of fracture of the wrist or after early removal of the plaster

Item	REF.OFT500 DX • REF.OFT500 SX
Size	universal
Length cm	15,5
Colour	blue

Specify right or left



Wrist immobiliser

- Aluminum thumb **SPLINT** preformed and pliable (removable)
- Aluminum palm **SPLINT** preformed and pliable (removable)
- Reinforcing dorsal plastic **STAYS**
- Double-lined foam **PADDING**
- **FASTENING** with Velcro® straps and buckles

INDICATIONS

- Wrist sprain
- Wrist tendonitis
- Consequences of wrist fracture or after early removal of the plaster
- Damage to ulnar section of metacarpophalangea of the thumb

Item	REF.OFT502 DX • REF.OFT502 SX
Size	universal
Length cm	16
Colour	blue

Specify right or left



Knee immobiliser

- **STRUCTURE** made of double-lined soft foam
- Internal **LINING** in terrycloth
- Rigid but moldable medial, side and back **RODS** made of light metal alloy
- Patella **OPENING**
- Adjustable supporting **PLATES**, inside and lateral

INDICATIONS

- Immobilisation with knee in extension position (post-surgical or post-trauma care)

Item	REF.OFT700/2	REF.OFT700/3
Immobilisation at	0°	0°
Length cm	20	30
Thigh circumference cm	24/35	33/45
Colour	blue	

Fits right and left



Bivalve ankle brace

- **EASY TO WEAR** thanks to full opening
- Velcro® adjustment **FITS ANY HEEL** width
- Double-lined foam **PADDING**
- Anatomically shaped lightweight durable plastic bivalve **STRUCTURE**

INDICATIONS

- Grade II and III sprains
- Functional recovery phase after surgical treatment
- Functional recovery after malleolar fractures
- Chronic instability of TT and sporting activity

Item	REF.OFT800
Size	universal
Height cm	16
Colour	black

Fits right and left



Walker

- **STRUCTURE** made from lightweight, durable plastic
- **STRUTS** held in place with Velcro® for a firm grip on the padding
- Velcro® fastening **PADDING**
- Inner **PADDING** of the sole is made from non-toxic Evazote®
- Non-slip rubber **SOLE**

INDICATIONS

- Conservative treatment of malleolar fractures
- Severe sprains grade 3
- Post-surgical phase after osteosynthesis of ankle fractures, stabilisation of ligament injuries of the ankle and Achilles tendon
- Stable foot, ankle, tibia and distal fibula fractures

Item	REF.OFT805		
Size	M	L	XL
Shoe size	23/25	26/28	29/31
Height from ground cm	22	25	27
Colour	black		

Fits right and left



Shower for lower leg and foot

- **FASTENS** with a Velcro® strap around the calf
- Polypropylene **SPLINT** can be heated and moulded to fit

INDICATIONS

- Foot drop
- Mild to moderate spastic equinus



REF.9391 R
REF.9391 L
semifinished of
REF.OFT810
Specify right or left

Item	REF.OFT810 • REF.9391 R • REF.9391 L				
Size	1.	2.	3.	4.	5.
Calf circumference cm	15/19	19/23	22/24	23/27	25/29
Foot section length cm	10	11	12	13	15
Splint height cm	16	19	23	25	28
Colour	semitransparent				

Specify right or left

ORTHO 14

Supporting the next generation



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**

Take care feel better