

Why choose flex-B™

1. Supports the spine
2. Ensures the right compression by increasing intra-abdominal pressure
3. Guarantees effectiveness, fit, and comfort thanks to the combined work of the elastic fabric and the metal frame
4. Reinforced and elastic, it is worn like a backpack


90CDP0011TEN02.1224

flex-B™

Padded dorso-lumbar orthopaedic orthosis with a reinforced structure
REF.33009

Code	REF.33009					
Size	XS	S	M	L	XL	XXL
Pelvic circumf. cm	70/80	80/90	90/100	100/110	110/120	120/130
Height at front cm	16	16	16	16	16	16
Back height cm	50	50	50	52	52	52
Colour	grey					

OPTION: Frame REF.33009/T size 1-L, this special frame fits sizes XS to M - with XXL back height



Indications

- Back pain caused by arthrosis or kyphosis
- Recovery from fractures of the thoracic and/or upper lumbar vertebrae
- Support for thoracic vertebral compression fracture caused by osteoporosis or metastasis
- Other conditions as prescribed by the doctor



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



flex-B™


Padded dorso-lumbar orthopaedic orthosis with a reinforced structure




**ORTHOSERVICE
RO+TEN**
Take care feel better



**ORTHOSERVICE
RO+TEN**
Take care feel better



flex-B™



**CARBONTRACE
FABRIC
STRUCTURE**
with static-free and
bacteriostatic **carbon fiber**



**REINFORCEMENT
FRAME IN
ALUMINUM**
lightweight, durable, and
easily moldable by the
orthopedic technician



**3+3 ABDOMINAL
STRAPS FOR
PERSONALIZED
SUPPORT**
reinforces
containment



**FRONT CLOSURE
WITH A REINFORCED
ABDOMINAL AREA**
equipped with a **thumbhole loop**

**ADJUSTABLE
AND ANTI-REMOVAL
STRAPS**
worn like a backpack



**REAR
ELASTIC
STRAP**
positionable to
position pressure
where needed

Easy to wear,
easy to carry.