

polFIT 22

Long-size immobilizer

with thumb support,
open



**ORTHOSERVICE
RO+TEN**

Take care feel better

polFIT 22

Long-size immobilizer with thumb support, open

Size

Colore

fits right and left

REF. PS502

universal

black

Size

Wrist circumference cm

Length cm

Colour

specify right or left

REF. P1121 specify right or left

S

13/16

M

16/18

L

18/20

XL

20/23

26

26

26

26

black



RADIO
LUCENT



EASY
FASTENING



3D FABRIC



BREATHABLE

Indications

- A compound fracture of the first metacarpal
- A sprain of the radial ulnar collateral ligament of the metacarpophalangeal joint
- A suspected scaphoid fracture

Contraindications

- Currently no known

Characteristics and materials

Polfit 21

- Full dorsal opening
- 3D breathable fabric for structure
- Four rings and Velcro® straps lock
- Dorsal reinforcing stays
- Aluminium rigid palmar stay, preformed and removable

PS502

- Cotton-blend fabric support with padding around the splints
- Adjustable Velcro® thumb strap
- Velcro® wrist strap



www.orthoservice.com



Headquarter: ORTHOSERVICE AG

Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89

info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH

Flugstraße 8 · D-76532 Baden-Baden · Deutschland

Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

Orthoservice Polska sp. z o. o.

ul. Warszawska 416a · 42-209 Częstochowa · Polska

Tel. +48 (0) 34 340 13 10

info@orthoservice.pl · www.orthoservice.pl

Sede italiana: RO+TEN s.r.l.

Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia

Sede operativa e amministrativa:

Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia

Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34

info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):

Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



**ORTHOSERVICE
RO+TEN**

Take care feel better