

# carbon afo

## Carbon-fibre ankle foot orthosis



**A lot of new features**  
for an even better result

### THE SPIRAL

has been reinforced  
to increase foot **support**

### THINNER CALF PADS

to improve the profile  
and fit the brace better



### STRENGTHENED ANTERIOR PLANTAR

for better  
toe support



### SUPERIOR COMFORT

due to the **softened joint**  
between the spiral  
and the sole

### FUNCTIONAL CALF GRIP

right mix of **flexibility**  
(to permit flexible  
adjustment  
to different calf  
circumferences) and  
**stiffness** (to avoid  
twisting during use)



**ORTHOSERVICE  
RO+TEN**

Take care feel better

# carbon afo

## Carbon-fibre ankle foot orthosis

REF. KTAFO

### Indications

- Foot drop
- Ictus
- Multiple sclerosis
- Neurological disorders
- Amputation of forefoot

### Contraindications

- Macroscopic deformity of the foot and ankle
- Acute oedema
- Moderate to severe spasticity
- Contracture of front and rear tibial muscles that cannot be healed by a simple positional device



RADIO  
LUCENT



LOW  
PROFILE



CARBON  
FIBER



SINGLE PATIENT  
MULTIPLE USE

Code

Size

Height cm

Shoe size

Colour

REF. KTAFO

S	M	L
36	39	42
33/38	38/43	43/48
black		

specify right or left



orthoservice.com



**Headquarter: ORTHOSERVICE AG**

Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89

info@orthoservice.com · www.orthoservice.com

**Niederlassung Deutschland: Orthoservice Deutschland GmbH**

Flugstraße 8 · D-76532 Baden-Baden · Deutschland

Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

**Siedziba w Polsce: Orthoservice Polska Sp. z o. o.**

ul. Warszawska 416a · 42-209 Częstochowa · Polska

Tel. +48 (0) 34 340 13 10

info@orthoservice.pl · www.orthoservice.pl

**Sede italiana: RO+TEN S.r.l.**

**Sede legale:** Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia

**Sede operativa e amministrativa:**

Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia

Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34

info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):

Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



**ORTHOSERVICE  
RO+TEN**

Take care feel better