

rizostab

Thumb spica



ORTHOSERVICE
RO+TEN

Take care feel better

rizostab

Thumb spica

REF. 21006



It is shorter than the thumb so it is easier for the thumb to grip with the forefinger



The aluminium core can be moulded without heat to fit the patient's hand



The Velcro® strap around the thumb can be adjusted

Completely covered with comfortable fabric

Sizing is easy to manage

Code	REF. 21006	
Size	S	L
Wrist circumference cm	14-18	18-22
Colour	grey	

Fits right and left



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com

What does it do?

It immobilises the trapeziometacarpal joint to relieve pain.

Indications

- Rhizarthrosis
- Acute or chronic joint pain
- Skier's thumb
- Tendonitis



LATEX FREE



RADIO LUCENT



CUSTOMIZABLE FIT



**ORTHOSERVICE
RO+TEN**

Take care feel better