

rizostab

**Thumb
spica**



**ORTHOSERVICE
RO+TEN**

Take care feel better

rizostab

Thumb spica

REF. 21006



It is shorter than the thumb so it is easier for the thumb to grip with the forefinger



The aluminium core can be moulded without heat to fit the patient's hand



The Velcro® strap around the thumb can be adjusted

Completely covered with comfortable fabric

Sizing is easy to manage



LATEX FREE



RADIO LUCENT



CUSTOMIZABLE FIT

Code

REF. 21006

Size

S

L

Wrist circumference cm

14-18

18-22

Colour

grey

Fits right and left

Indications

- Rhizarthrosis
- Acute or chronic joint pain
- Skier's thumb
- Tendonitis

What does it do?

It immobilises the trapeziometacarpal joint to relieve pain.

ORTHOSERVICE
RO+TEN

Eco
attentive



UNI EN ISO 13485:2016
UNI EN ISO 9001:2015

Orthoservice AG
www.orthoservice.com



Headquarter: ORTHOSERVICE AG

Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland

Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89

info@orthoservice.com · www.orthoservice.com

Niederlassung Deutschland: Orthoservice Deutschland GmbH

Flugstraße 8 · D-76532 Baden-Baden · Deutschland

Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13

info@orthoservice.de · www.orthoservice.de

Orthoservice Polska sp. z o. o.

ul. Warszawska 416a · 42-209 Częstochowa · Polska

Tel. +48 (0) 34 340 13 10

info@orthoservice.pl · www.orthoservice.pl

Sede italiana: RO+TEN s.r.l.

Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia

Sede operativa e amministrativa:

Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia

Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34

info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):

Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland



ORTHOSERVICE
RO+TEN

Take care feel better



orthoservice.com