

## Several applications, several benefits

APPLICATIONS	BENEFITS
post-arthroscopy care	oedema prevention, cryotherapy and movement control
conservative treatment of hip dislocations	movement control (flexion and extension and especially postoperative abduction/adduction)
post-surgical care (also for replacements or revision surgery)	
arthrosis of the hip joint	pain relief
coxalgia or inoperable joint misalignment	

INNOVATION 2

Only 3 fits for any request

The versatility of the brace reduces the need to have to keep stocks in stock, knowing nonetheless that every client's request will be satisfied.

Size	S	M	L
Pelvis circ. cm	72/104	104/137	137/170
Proximal thigh circ. cm	40/58	54/74	68/90

9000P0155EN001.0322

## hipocross 2.0

### Brace for hip and coxalgia

REF.91006N - REF.91007N

#### Indications

##### REF.91006N:

- Post-arthroscopy care.
- Post-surgical care (compatible with the degree of stabilization required).
- Conservative treatment of hip dislocations.
- Arthrosis of the hip joint.
- Coxalgia or inoperable joint misalignment.



Item	REF.91006N (Hipocross 2.0)		
Size	S	M	L
Pelvis circ. cm	72/104	104/137	137/170
Proximal thigh circ. cm	40/58	54/74	68/90
Colour	black		

specify right or left

Item	REF.91007N (Adjustable telescopic rod for Hipocross 2.0)		
Size	universal		
Colour	black		

specify right or left



orthoservice.com



Headquarter: ORTHOSERVICE AG  
Via Milano 7 - CH-6830 Chiasso (TI) - Switzerland  
Tel. 0041 (0) 91 822 00 88 - Fax 0041 (0) 91 822 00 89  
info@orthoservice.com - www.orthoservice.com  
Niederlassung Deutschland: Orthoservice Deutschland GmbH  
Flugstraße 8 - D-76532 Baden-Baden - Deutschland  
Tel. 0049 (0) 7221 991 39 11 - Fax 0049 (0) 7221 991 39 13  
info@orthoservice.de - www.orthoservice.de  
Siedziba w Polsce: Orthoservice Polska Sp. z o. o.  
ul. Warszawska 416a - 42-209 Częstochowa - Polska  
Tel. +48 (0) 34 340 13 10  
info@orthoservice.pl - www.orthoservice.pl  
Sede italiana: RO+TEN S.r.l.  
Sede legale: Via Marco De Marchi, 7 - I-20121 Milano (MI) - Italia  
Sede operativa e amministrativa:  
Via Comasina, 111 - I-20843 Verano Brianza (MB) - Italia  
Tel. 0039 039 601 40 94 - Fax 0039 039 601 42 34  
info@roplusten.com - www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):  
Orthoservice AG (CH) - 6830 Chiasso (TI) - Switzerland



REF.91007N  
Adjustable (telescopic) rod for Hipocross



hipocross 2.0

## Brace for hip and coxalgia



ORTHOSERVICE  
RO+TEN  
Take care feel better



ORTHOSERVICE  
RO+TEN  
Take care feel better

# Already a winner, but even better from today: 3 absolute innovations

## TELESCOPIC AND ARTICULATED ROD

Allows numerous regulations:

- **Height**, by sliding the two segments to the desired measurement.
- **Adduction**, thanks to the joint present at trochanteric height
- **Flexion-extension**, through a special graduated joint (from 0° to 120°).

INNOVATION 1



INNOVATION 3

## VELCRO® FASTENING SYSTEM FOR THE STRAPS ON BOTH THE THIGH AND THE HIP BONE

The Velcro® allows quick and millimeter precision fastening.



## ANTI-EDEMA CONTAINMENT FABRIC AND CRYOTHERAPY GEL

The compression exerted by the combination of the special containment fabric and the gel (that has to be inserted in the pocket in the area of the hip), help the reabsorption of the edema. The gel is cooled in the freezer.



## RIGID BAR

A rigid bar that locks the joint in flexion and neutral abduction during the initial days following a dislocation.



## X-SHAPED TIE RODS ON THE SIDE FOR EASY ADDUCTION

The 2 robust rods that cross on the side allow the dynamic adjustment of the adduction.



Attention to detail makes the difference