

## Reasons for choosing this brace

- 1 **Practicality, stability, and maximum comfort in a low-profile orthosis**
- 2 **Lanyard and pulley traction system provides simple closure with minimum effort**
- 3 **Two sizes only – small and big- for a wide range of measurements**

4 **Modular, for different pathologies**



Lumbofast Pad  
Additional pads kit



9000P00179EN001723

## lumbofast

REF. 23070 Orthosis with a facilitated traction system to stabilize the lumbar rachis and relieve pressure with rigid thermoformable pads

Item	REF.23070	
Size	1	2
Pelvic circumference cm	70/106	106/145
Front height cm	16	16
Back height cm	24	24
Color	grey	

Option: REF.23070P additional pads kit for Lumbofast



### Indications

#### LUMBOFAST REF. 23070

Involved spine section L2-L5

- Lumbago
- Sciatic lumbalgia and cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas of the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- After-effects of stable compression fractures of the lumbar vertebrae
- Treatment of fracture of transverse processes

#### LUMBOFAST REF.23070 + LUMBOFAST PAD REF.23070P

Involved spine section L1-L5

- Lumbago
- Sciatic lumbalgia and cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas of the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- After-effects of stable compression fractures of the lumbar vertebrae
- Treatment of fracture of transverse processes
- Degenerative lumbar scoliosis
- Spondylosis with or without vertebral listhesis



Lumbofast

lumbofast

**Orthosis with a facilitated traction system**  
to stabilize the lumbar rachis and relieve pressure



orthoservice.com



Headquarter: ORTHOSERVICE AG  
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland  
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89  
info@orthoservice.com · www.orthoservice.com

ORTHOSERVICE DEUTSCHLAND GmbH  
Flugstraße 8 · D-76532 Baden-Baden · Deutschland  
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13  
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.  
ul. Warszawska 416a · 42-209 Częstochowa · Polska  
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785  
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.  
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia  
Sede operativa e amministrativa:  
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia  
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34  
info@roplusten.com · www.roplusten.com



ORTHOSERVICE  
RO+TEN  
Take care feel better



ORTHOSERVICE  
RO+TEN  
Take care feel better

**POLYPROPYLENE REINFORCEMENT**  
to protect the lanyard

**WIDER FRONT LEG OPENING**  
for an optimal fit

**LIGHTWEIGHT WITH A LOW PROFILE**

**ANATOMICAL DESIGN**



**INNOVATION**

**LANYARD AND PULLEY TRACTION SYSTEM**  
multiply the pulling force by 4 adapting the brace to the patient's conformation



**LANYARD TERMINALS**  
have an anatomical and easy grip



**FRONT HANDLES**  
for easy Velcro® closure



**BI-COUPLED BASE FABRIC**  
and breathable

**INNOVATION**

**ONLY TWO SIZES**  
to dress patients with pelvic circumference from 70 cm to 145 cm reducing stocks



The cross-use orthosis, even when it comes to heavy physical work