

lumbofast

Orthosis with a facilitated traction system

to stabilize the lumbar rachis and relieve pressure



**ORTHOSERVICE
RO+TEN**


Take care feel better

Reasons for choosing this brace

- 1 **Practicality, stability, and maximum comfort in a low-profile orthosis**
- 2 **Lanyard and pulley traction system provides simple closure with minimum effort**
- 3 **Two sizes only – small and big- for a wide range of measurements**

4

Modular, for different pathologies




Involved spine section L2-L5



Lumbofast

Lumbofast Pad
Additional pads kit



Involved spine section L1-L5



lumbofast

**POLYPROPYLENE
REINFORCEMENT**

to protect the lanyard

**WIDER FRONT
LEG OPENING**

for an optimal fit

**LIGHTWEIGHT
WITH A LOW PROFILE**

**ANATOMICAL
DESIGN**



INNOVATION

LANYARD AND PULLEY TRACTION SYSTEM

multiply the pulling force by 4 adapting the brace to the patient's conformation



LANYARD TERMINALS

have an anatomical and easy grip



FRONT HANDLES

for easy Velcro® closure

The cross-use
orthosis, even when
it comes to heavy
physical work



BI-COUPLED BASE FABRIC
and breathable

INNOVATION

ONLY TWO SIZES

to dress patients with pelvic circumference
from 70 cm to 145 cm reducing stocks



lumbofast

REF. 23070 Orthosis with a facilitated traction system to stabilize the lumbar rachis and relieve pressure with rigid thermoformable pads

Item	REF.23070	
Size	1	2
Pelvic circumference cm	70/106	106/145
Front height cm	16	16
Back height cm	24	24
Color	grey	

Option: REF.23070P additional pads kit for Lumbofast



ANATOMICAL
FIT



BREATHABLE



3D FABRIC



REGISTERED
EU, CH
AND UK DESIGN



SINGLE PATIENT
MULTIPLE USE



RIGID
THERMOFORMABLE
PAD



EASY
PULLING



Lumbofast

Indications

LUMBOFAST REF. 23070

Involved spine section L2-L5

- Lumbago
- Sciatic lumbalgia and cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas of the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- After-effects of stable compression fractures of the lumbar vertebrae
- Treatment of fracture of transverse processes

LUMBOFAST REF.23070 + LUMBOFAST PAD REF.23070P

Involved spine section L1-L5

- Lumbago
- Sciatic lumbalgia and cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas of the lumbosacral column and paravertebral lumbar contractures
- Postoperative recovery
- After-effects of stable compression fractures of the lumbar vertebrae
- Treatment of fracture of transverse processes
- Degenerative lumbar scoliosis
- Spondylosis with or without vertebral listhesis



orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**

Take care feel better