

easyco

Ribbed elastic orthosis



ORTHOSERVICE
RO+TEN

Take care feel better

easyco

Why chose Easyco?

- Easy to fit
- Easy to wear
- 2 **Basic** and **Plus** versions (with tie rods) for every need

EASYCO BASIC

FRONT
STAYS

FRONT
HANDLE

to facilitate closing

FRONT
VELCRO®
CLOSURE

Say hello Easyco
and goodbye to
back pain!

EASYCO PLUS

**CRISS-CROSSED
EXTERNAL
ELASTIC
TIE-RODS**

as a reinforcement,
with a Velcro®
fastening system
(only plus version)



**REINFORCED
ELASTIC BACK**

(only basic version)



**4 LIGHTWEIGHT AND
FLEXIBLE RODS ON THE BACK**

(width 20 mm, height 20/22 cm)

**STRUCTURE MADE OF
RIBBED ELASTOMER
FABRIC**



**2 DIFFERENT COLORS:
BEIGE AND BLACK**

easyCO basic

Ribbed elastic orthosis

REF. 23770B (beige) – REF. 23770S (black)

Indications

- Paravertebral muscle contractures
- Mild traumas
- Mild lumbago

easyCO plus

Ribbed elastic orthosis with tie-rods

REF. 23870B (beige) – REF. 23780S (black)

Indications

- Paravertebral muscle contractures
- Mild traumas
- Mild lumbago
- Mild sciatic lumbalgia and cruralgia
- Post-operative recovery



SINGLE PATIENT
MULTIPLE USE



BREATHABLE



Easyco Basic



Easyco Plus

Item	REF.23770B (beige) • REF.23770S (black) / REF.23870B (beige) • REF.23870S (black)					
Size	XS	S	M	L	XL	XXL
Waist circumference cm	50-60	60-75	75-90	90-105	105-120	120-135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26
Colour	black/beige					



orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**

Take care feel better