



When to propose the Stepper Line

The Stepper Orthoses have been designed to accommodate a **swollen and painful foot after surgery**. They provide support where needed while also offering **comfort and reassurance**. They protect injuries from accidental impacts, lighten the load, and instill a **pleasant sense of security** for walking stability.

910CP0102021INT02.0226

stepper

M23001 • Stepper SURG post-op orthosis for the foot

Code	REF.23001						
Size	3XS	XXS	XS	S	M	L	XL
Shoe size	25/32	33/36,5	37/38,5	39/41	41,5/43	43,5/45	45,5/47

Supplied individually- ambidextrous

Indications

- Post-op for stabilization and reduce load from the foot (e.g. foot osteotomies, hallux valgus,...)

M23002 • Stepper DIGIT post-op orthosis to reduce load on the forefoot

Code	M23002					
Size	XS	S	M	L	XL	
Shoe size	33,5/37	37,5/39,5	40/41,5	42/44	44,5/47	

Supplied individually- ambidextrous

Indications

- Post-op orthosis to reduce load on the forefoot

M23003 • Stepper CAL Post-op orthosis to reduce load on the hindfoot

Code	M23003					
Size	XS	S	M	L	XL	
Shoe size	33,5/37	37,5/39,5	40/41,5	42/44	44,5/47	

Supplied individually- ambidextrous

Indications

- Post-operative recovery after both soft tissue and bony heel surgery
- Selective pressure relief after fractures
- Restoration of hindfoot ulcers

M23004 • Toe cap for Stepper SURG and Stepper DIGIT

Code	M23004					
Size	XS	S	M	L	XL	

Supplied individually- ambidextrous



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. +49 (0) 91 822 00 88 · Fax +49 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. +49 (0) 7221 40 533 0 · Fax +49 (0) 7221 40 533 20
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



Stepper SURG



Stepper DIGIT



Stepper CAL



Stepper toe cap



RADIO LUCENT

RIGID SOLES

SINGLE PATIENT MULTIPLE USE

UNISEX

stepper

Post-op orthoses for foot



ORTHOSERVICE
RO+TEN

Take care feel better



ORTHOSERVICE
RO+TEN

Take care feel better

Why choose these orthoses for post-op

- 1 **Maximum functionality** thanks to the ultra-rigid sole
- 2 **Maximum protection** thanks to the squared-off front shape to protect the toes
- 3 **Maximum comfort** thanks to the lightweight and breathable upper and the wide opening for fit
- 4 **Unisex**, which the warehouse thanks for



UNISEX



stepper SURG

- No height compensation necessary on the contralateral limb, low-profile
- Reduction of pressure on the heel and forefoot thanks to the rigid sole
- Boat-shaped profile that facilitates foot roll-off

stepper toe cap

- Toe cap for Stepper SURG and Stepper DIGIT



stepper DIGIT

- Wedge at 15° to reduce load on the forefoot



stepper CAL

- Selective unloading of the pressure on the heel



FULL EVA INSOLE

TOTAL OPENING

UPPER

- Lightweight and breathable
- Rear counter for better hindfoot stabilization
- Fully openable for easy fitting
- Velcro® strap on the instep with reinforcement ring; with protective padding
- Symmetric closure on the forefoot, without buckles

RIGID SOLE

- Stiff non-slip sole with reinforcement bar
- Square-shaped front for toe protection
- Removable EVA insole to be optionally adapted or replaced

Exclusive design by Orthoservice RO+Ten