

Walker with two inflatable air padding



airstep tight walker

REF.550N



airstep tight walker
(short)

REF.560N

INDICATIONS

- Conservative and post-operative treatment of foot/ankle and leg fractures
- Post-operative immobilisation after surgical treatment of foot/ankle fractures for patients with short limbs (Airstep tight Walker short)
- Serious bone distortion grade 1 to 3
- Conservative treatment after a tear in the Achilles tendon
- Post-surgical treatment after tenorrhaphy of the Achilles tendon
- Pre- and post-operative/post-traumatic immobilisation

CONTRAINDICATIONS

The application of this brace is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments

SELECTION/SIZE

Item	REF.550N/560N(short)			
Size	S	M	L	XL
Shoe size	35/39	39/44	44/47	>47
Colour	grey			

fits right and left



AIR PAD



RADIO LUCENT



SINGLE PATIENT MULTIPLE USE



www.orthoservice.com



www.kiwa.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**
Take care feel better

airstep tight walker

Walker with two
inflatable air padding

NEW 2023



**ORTHOSERVICE
RO+TEN**
Take care feel better

The new and improved Walker



normal version



short version

- 1 Reinforcing shells** on the front and back, for increased stability
- 2 Inflatable air pads** for optimum fit to the limb and improved stability
- 3 Structure in lightweight, sturdy plastic**; terminal part in rubbery material
- 4 Integrated pump** to inflate/deflate the padding
- 5 Sole** in non-slip material
- 6 Optional padding** for additional adaptation to the leg

Can be put on in just five steps



OPTION:



REF.505
Diabetic Set, fits right and left.



REF.501
Wedges for walker specify right and left.

Two versions for every need