

epifeel

REF.24022

Epicondylitis and epitrochleitis elbow strap with pressure pad

Code	REF.24022
Size	universal
Colour	grey

Fits right and left

OPTION: REF.ICEGEL Pad for cold therapy with Velcro® fastening system

EASY
PULLINGRADIO
LUCENTLOW
PROFILESINGLE PATIENT
MULTIPLE USE

Indications

- Acute or chronic epicondylitis and epitrochleitis of the elbow
- Muscle contractures
- Prevention of inflammatory processes during sports activity

For epicondylitis

For epitrochleitis

epifeel

Epicondylitis and epitrochleitis elbow strap with pressure pad



orthoservice.com


UNIVERSITÀ DI TRIESTE
SISTEMA DI GESTIONE
CERTIFICATO
UNI EN ISO 9001:2015

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



**ORTHOSSERVICE
RO+TEN**
Take care feel better



**ORTHOSSERVICE
RO+TEN**
Take care feel better

Why choose Epifeel

- Designed for those suffering from epicondylitis or epitrochleitis
- Comfortable and invisible under clothes
- Lightweight and easy to adjust

ELASTIC ELEMENT AT THE LEVEL OF THE CLOSURE LOOP

to avoid the lace effect

LYCRA PADDING COUPLED



PRESSURE PAD TO FOCUS THE THRUST

made of non-slip TPE for an easy positioning of the brace

UNIVERSAL SIZE

POLYAMIDE STRUCTURE

for optimal pressure distribution

Ideal for sports and daily activities

VELCRO® FASTENER



ADJUSTABLE STRAP

for accurate compression control



ADJUSTMENT OF CORRECT TENSION

thanks to the white reference inside the ring

When propose Epifeel

From tennis elbow to mouse posture; **Epifeel** is used to relieve pain by unloading the inflamed tendon.