

Why choose this ankle brace

- 1 Easy to wear thanks to open heel closure
- 2 Simple and quick drawstring adjustment system
- 3 Anatomical and low-profile design for use with most shoes
- 4 Multifunctional use, even during sports activities

When suggest Malleorapid

Malleorapid is recommended for **II degree sprain trauma** and as an adjuvant in the resumption of **sporting activity**.

The combined action of the steel stays and the 8-figure elastic straps limits inversion/eversion of the ankle

9000P00193EN001.0124

malleorapid

REF. 24081
Ankle brace with steel stays and figure-8 elastic straps

Code	REF. 24081		
Size	1	2	3
Circumference ankle cm	19-22	22-26	26-30
Colour	grey		

Fits right and left



Indications

- Conservative treatment following a grade 2 sprain trauma
- Acute tendinopathies
- Post-surgical rehabilitation
- Aid to resuming sporting activity after grade 2 and 3 sprain trauma
- Chronic joint pain caused by tibiotarsal arthrosis



orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

ORTHO SERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHO SERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



malleorapid

Ankle brace with steel stays and figure-8 elastic straps



ORTHO SERVICE RO+TEN
Take care feel better

ORTHO SERVICE RO+TEN
Take care feel better

ANATOMIC DESIGN

to perfectly adapt to the morphology of the ankle

MOLDABLE MEDIAL AND LATERAL STAYS MADE FROM A ROUND STEEL ROD

to increase the ankle stability



ADDITIONAL LATERAL AND MEDIAL REINFORCEMENT FABRIC

to protect against contact with the stays



DRAWSTRING ADJUSTMENT SYSTEM WITH QUICK-RELEASE BUCKLE

for fast and precise adaptation



TIGHTNESS AND MAXIMUM BREATHABILITY

thanks to the three-dimensional mesh fabric

LYCRA POCKET

for storing the lanyard



OPEN HEEL WITH VELCRO® CLOSURE

to put on and take off quickly the brace



ELASTIC CROSSOVER STRAPS

with Velcro® closure on the instep for maximum stability



The ankle brace for training safely or for daily activities