

Why choose this ankle brace

- Easy to wear thanks to open heel closure
- Simple and quick drawstring adjustment system
- Anatomical and lowprofile design for use with most shoes
- **Multifunctional use**, even during sports activities

When suggest **Malleorapid**

Malleorapid is recommended for **II degree sprain trauma** and as an adjuvant in the resumption of **sporting activity**.

The combined action of the steel stays and the 8-figure elastic straps limits inversion/eversion of the ankle





malleorapid

Ankle brace with steel stays and figure-8 elastic straps

Code	REF. 24081			
Size	0	1	2	3
Circumference ankle cm	16-19	19-22	22-26	26-30
Colour	grey			

Fits right and left











BREATHABLE

- Conservative treatment following a grade 2 sprain trauma
- Acute tendinopathies
- Post-surgical rehabilitation
- Aid to resuming sporting activity after grade 2 and 3 sprain trauma
- Chronic joint pain caused by tibiotarsal arthrosis















Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland Tel. 0041 (0) 91 822 00 88 Fax 0041 (0) 91 822 00 89 info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH

Flugstraße 8 · D-76532 Baden-Baden · Deutschland Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13 info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.

ul. Warszawska 416a · 42-209 Częstochowa · Polska Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785 info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd

Units 47&50 North East BIC Wearfield Sunderland Enterprise Park East Sunderland SR5 2TA Tel +44 (0) 191 51 66 220

info@orthoservice.uk · www.orthoservice.uk

Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia Sede operativa e amministrativa:

Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia Tel. 0039 0362 18 51 200

info@roplusten.com · www.roplusten.com



