


elbo PRO

Telescopic
elbow brace
with graduated
movement



ORTHOSERVICE
RO+TEN

Take care feel better



...satisfies all sizes in case of:

- **Compound fractures of the elbow, distal humerus, proximal radius and ulna**
- **post-operative immobilisation following tendon and ligament repairs**
- **following surgical and conservative treatment of elbow fractures**
- **elbow subluxations and dislocations**
- **chronic elbow injuries**
- **elbow hyperextension**
- **post-traumatic treatment**

Efficient and versatile...

TELESCOPIC

Easily adapts to every arm eliminating
the need for different sizes

THE JOINT IS EASY AND SAFE

It is regulated without the need for tools
but cannot be unintentionally released




TELESCOPIC LENGTH ADJUSTMENT OF THE ARTICULATED HINGE

thanks to easy-to-use buttons
on the valves



FIRM ANCHORAGE

(the brace doesn't rotate) thanks
to the **polyamide valves** with
mouldable aluminium **arm grip**



EASY FLEXION-EXTENSION ADJUSTMENT

thanks to the **'Pull&Set'** joint with spring-loaded sliders (pull outwards, position to desired degree, release)



MAXIMUM COMFORT

thanks to the **heat-sealed padding** in soft lycra that comes into contact with the skin

Elbo PRO,
the new era
for elbow braces

Included accessories

LENGTH-ADJUSTABLE SHOULDER STRAP

with **padded shoulder protection** to be worn on the ipsilateral shoulder for optimal comfort



THE SET FLEXION-EXTENSION

cannot be altered by arm movements or impacts against surfaces

NEUTRAL HAND SUPPORT

for forearm rotation control and a feeling of stability

elbo PRO

M 24017R

M 24017L

Telescopic elbow brace with graduated movement

Code	M 24017R • M 24017L
Size	universal
Length of the hinges cm	36-46
Colour	dark grey

Right or left

CUSTOMIZABLE
FITEASY
ADJUSTING

TELESCOPIC

SINGLE PATIENT
MULTIPLE USE

Indications

- Compound fractures of elbow, distal humerus, proximal radius and ulna
- Post-operative immobilisation following tendon and ligament repairs
- Useful in the rehabilitation phase following surgical and conservative treatment of elbow fractures
- Elbow subluxations and dislocations
- Chronic elbow injuries
- Elbow hyperextension
- Post-traumatic treatment



orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com



**ORTHOSERVICE
RO+TEN**

Take care feel better