

hernix sanitary underwear


Hernix line
for intimate relief



ORTHOSERVICE
RO+TEN

Take care feel better

The full range of hernia support and containment briefs



The most amazing effect of the Hernix **hernia support and containment line** of briefs, is the fact of always feeling protected and in good shape. Thanks to **meticulous design, exclusive fabrics**, and the choice of **popular colors**, the Hernix briefs are discrete and suitable **for any activity**, without the fear of being noticed. They sustain and compress where needed while letting the skin breathe. That is why they **can be worn all day long**.

hernix sanitary underwear

ITEM	TYPE		CONTAINMENT			HEIGHT			COLOR	INCLUDED
	BRIEF	BOXER	STRONG	MEDIUM-STRONG	LIGHT	HIGH	MEDIUM	LOW		
HERNIA BRIEF										
Hernix 10	●		●			●			light grey	M24699M
Hernix 42	●			●			●		anthracite grey	M24699M
Hernix 43A	●			●			●		anthracite/ light grey	M24699M
Hernix 43C	●			●			●		light grey/ anthracite	M24699M
Hernix 49		●		●			●		anthracite grey	M24699M
Hernix 90	●		●			●			nude	M24699W
CONTAINMENT BRIEFS										
Hernix 45	●			●			●		anthracite grey	
Hernix 46	●				●		●		anthracite grey	
Hernix 47					●			●	anthracite grey	

HERNIA BRIEFS

What are they for

They are used to contain hernias where they do not exist contraindications.

What are they like

Thin and soft against the skin, they are made of exclusive cotton-based fabrics. Thanks to the appropriately shaped reinforcements, they guarantee effective containment of the hernia. Beautiful to look at and pleasant to wear.

How do they make you feel

They provide an indispensable feeling of security, freedom of movement and comfort, in any situation.

CONTRAINDICATIONS

Do not use in case of incarcerated or strangulated hernias, at an advanced stage of inflammation or irreducible.

Code	M24610 • M24642 • M24643A • M24643C • M24649 • M24690							
Size	1	2	3	4	5	6	7	8
Hip circumf. cm	83-87	88-92	93-97	98-103	104-109	110-115	116-122	123-130

INCLUDED: M24699M HERNIX P Couple of pads for Hernix, men
M24699W HERNIX P Couple of pads for Hernix, female

hernix 10

Classic hernia brief



Sizes:
1-8



Color:
Light grey



Height:
High




Containment:
Strong



Type:
Brief



Composition:
82% CO
13% EA
5% PA



INCLUDED:
M24699M HERNIX P
Couple of pads for Hernix, men

M24610

INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

CHARACTERISTICS

- Very supportive cotton fabric structure
- Reinforced and shaped internal structure
- Internal under-thigh pullers
- Soft, stretchable front with side openings
- Fitted for hernia pads (M24699M)

Containment hernia brief



Sizes:
1-8



Color:
Anthracite grey



Height:
Medium



Containment:
Medium-strong



Type:
Brief



Composition:
**82% CO · 10% EA
6% PA · 2%PL**



Two-tone containment hernia brief



Sizes:
1-8



Color:
**Anthracite/ light grey (43 A)
Light grey/ anthracite (43 C)**



Height:
Medium




Containment:
Medium-strong



Type:
Brief



Composition:
**82% CO · 10% EA
6% PA · 2%PL**



INCLUDED:
M24699M HERNIX P

Couple of pads for Hernix, men


M24642

INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

CHARACTERISTICS

- Hypoallergenic, thin, soft and comfortable cotton fabric structure
- Reinforced and shaped inner structure
- Under-thigh inner pullers, with calibrated traction
- Closed and shaped front, lined with soft cotton jersey
- Fitted for hernia pads (M24699M)



INCLUDED:
M24699M HERNIX P

Couple of pads for Hernix, men

**M24643A
M24643C**

INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

CHARACTERISTICS







- Hypoallergenic, thin, soft and comfortable cotton fabric structure
- Reinforced and shaped inner structure
- Under-thigh inner pullers, with calibrated traction
- Closed and shaped front, lined with soft cotton jersey
- Fitted for hernia pads (M24699M)


Hernia boxer briefs

-  Sizes: **1-8**
-  Color: **Anthracite grey**
-  Height: **Medium**
-  Containment: **Medium-strong**
-  Type: **Boxer**
-  Composition: **68% CO • 17% PA
13% EA • 2% PL**



Female hernia brief

-  Sizes: **1-8**
-  Color: **Nude**
-  Height: **High**
-  Containment: **Strong**
-  Tipologia: **Female brief**
-  Composition: **82% CO • 13% EA
5% PA**



INCLUDED:
M24699M HERNIX P
Couple of pads for Hernix, men


M24649

INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

CHARACTERISTICS

- Hypoallergenic, thin, soft and comfortable cotton fabric structure adapted for boxer
- Reinforced and shaped inner structure
- Under-thigh inner pullers, with calibrated traction
- Closed and shaped front, lined with soft cotton jersey
- Fitted for hernia pads (M24699M)



INCLUDED:
M24699W HERNIX P
Couple of pads for Hernix, female

M24690

INDICATIONS

- Inguinal hernias, ptosis, and prolapse containment
- Post-op
- Hernia recurrence prevention

CHARACTERISTICS

- Very supportive, hypoallergenic, soft and comfortable cotton fabric structure.
- Reinforced and shaped internal structure for correct support of the abdominal wall
- Soft finishing on the leg area for maximum comfort
- Fitted for hernia pads (M24699W)
- Equipped with a gusset opening

CONTAINMENT BRIEFS

hernix 45

What is their purpose

Hernix 45: provide support and containment to the abdominal area following surgery or in cases of tissue relaxation.

Hernix 46 and Hernix 47: contain and support the abdominal and genital area during sporting activity or in any situation that causes severe strain.

How do they feel

They contain efficiently, are comfortable and breathable. You put them on and forget you are wearing them.

How do they make you feel

In good shape and well-protected, ready to perform in sports activities or any activity.

CONTRAINDICATIONS

Do not use in case of incarcerated or strangulated hernias, at an advanced stage of inflammation or irreducible.



Post op containment briefs



Size:
1-8



Color:
Anthracite grey



Height:
Medium



Containment:
Medium-strong



Type:
Brief



Composition:
**82% CO • 10% EA
6% PA • 2% PL**



M24645

INDICATIONS

- Post-abdominal surgery or in case of tissue relaxation for abdominal wall containment
- Genital support

CHARACTERISTICS

- Soft, thin, supportive cotton fabric structure
- Reinforced and shaped inner structure for effective containment of the abdominal wall and support of the genital region
- Closed, soft, highly stretchable waistband

Code	M24645 • M24646 • M24647							
Size	1	2	3	4	5	6	7	8
Hip circumf. cm	83-87	88-92	93-97	98-103	104-109	110-115	116-122	123-130

Sports containment briefs



Size:
1-8



Color:
Anthracite grey



Height:
Medium



Containment:
Light



Type:
Brief



Composition:
**82% CO • 10% EA
6% PA • 2% PL**



M24646

INDICATIONS

- Preventive action genital support during work activities and sports

CHARACTERISTICS

- Medium containment cotton fabric structure, thin and soft
- Closed front, with calibrated containment for correct containment of the genital region

Jockstrap



Size:
1-8



Color:
Anthracite grey



Height:
Low



Containment:
Light



Type:
Jockstrap



Composition:
**65% PA • 18% EA
15% CO • 2% PL**



M24647

INDICATIONS

- Preventive action genital support during work activities and sports

CHARACTERISTICS

- Tailored cotton fabric front for correct containment of the genital region
- Soft elastic under thigh
- Soft finish crotch

hernix

sanitary underwear



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
 Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
 Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOSERVICE DEUTSCHLAND GmbH
 Flugstraße 8 · D-76532 Baden-Baden · Deutschland
 Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de

ORTHOSERVICE POLSKA Sp. z o.o.
 ul. Warszawska 416a · 42-209 Częstochowa · Polska
 Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
 Units 47&50 North East BIC
 Wearfield
 Sunderland Enterprise Park East
 Sunderland SR5 2TA
 Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
 Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
 Sede operativa e amministrativa:
 Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
 Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



ORTHOSERVICE
RO+TEN

Take care feel better