

**hernix** compression underwear

**Hernix line**  
for intimate relief



**ORTHOSERVICE**  
**RO+TEN**

Take care feel better

# The Hernix range of support and containment briefs



**T**he meticulous design and use of exclusive fabrics allows the Hernix range of briefs to **provide support where needed**. The wearer will feel protected and confident **during any activity** as Hernix sustains and compresses while also allowing airflow to the skin. Discreet and breathable, Hernix **can be worn all day long with comfort**.

# hernix compression underwear

ITEM	TYPE		CONTAINMENT			HEIGHT			COLOUR	PADS
	BRIEF	BOXER	STRONG	MEDIUM-STRONG	LIGHT	HIGH	MEDIUM	LOW		
<b>HERNIA BRIEF</b>										
<b>Hernix 10</b>	●		●			●			light grey	<b>M24699M</b>
<b>Hernix 42</b>	●			●			●		anthracite grey	<b>M24699M</b>
<b>Hernix 43A</b>	●			●			●		anthracite/ light grey	<b>M24699M</b>
<b>Hernix 43C</b>	●			●			●		light grey/ anthracite	<b>M24699M</b>
<b>Hernix 49</b>		●		●			●		anthracite grey	<b>M24699M</b>
<b>Hernix 90</b>	●		●			●			nude	<b>M24699W</b>
<b>CONTAINMENT BRIEFS</b>										
<b>Hernix 45</b>	●			●			●		anthracite grey	
<b>Hernix 46</b>	●				●		●		anthracite grey	
<b>Hernix 47</b>					●			●	anthracite grey	

# HERNIA BRIEFS

## What does Hernix do?

They are used to contain hernias where there are no contraindications.

## How do they make you feel?

Comforting sense of security, freedom and relaxed in any situation.

## What are they like to wear?

Made of exclusive cotton-based fabric weaves they are soft against the skin, aesthetically pleasant and easy to wear.

## CONTRAINDICATIONS

Not to be used if the hernia is incarcerated, strangulated, irreducible or at an advanced stage of infection.

Code	M24610 • M24642 • M24643A • M24643C • M24649 • M24690							
Size	1	2	3	4	5	6	7	8
Hip circumf. cm	83-87	88-92	93-97	98-103	104-109	110-115	116-122	123-130

**INCLUDED:** M24699M HERNIX P Couple of pads for Hernix, men  
M24699W HERNIX P Couple of pads for Hernix, female

hernix 10

## Classic hernia brief



Sizes:  
**1-8**



Colour:  
**Light grey**



Height:  
**High**




Containment:  
**Strong**



Type:  
**Brief**



Composition:  
**82% CO**  
**13% EA**  
**5% PA**



**INCLUDED:**  
**M24699M HERNIX P**

Couple of pads for Hernix, men

**M24610**

## INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

## CHARACTERISTICS

- Very supportive cotton fabric structure
- Reinforced and shaped internal structure
- Internal under-thigh pullers
- Soft, stretchable front with side openings
- Fitted for hernia pads (M24699M)

## Containment hernia brief



Sizes:  
**1-8**



Colour:  
**Anthracite grey**



Height:  
**Medium**



Containment:  
**Medium-strong**



Type:  
**Brief**



Composition:  
**82% CO · 10% EA  
6% PA · 2%PL**



## Two-tone containment hernia brief



Sizes:  
**1-8**



Colour:  
**Anthracite/ light grey (43 A)  
Light grey/ anthracite (43 C)**



Height:  
**Medium**




Containment:  
**Medium-strong**



Type:  
**Brief**



Composition:  
**82% CO · 10% EA  
6% PA · 2%PL**



**INCLUDED:**  
**M24699M HERNIX P**

Couple of pads for Hernix, men


**M24642**

### INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

### CHARACTERISTICS

- Hypoallergenic, thin, soft and comfortable cotton fabric structure
- Reinforced and shaped inner structure
- Under-thigh inner pullers, with calibrated traction
- Closed and shaped front, lined with soft cotton jersey
- Fitted for hernia pads (M24699M)



**INCLUDED:**  
**M24699M HERNIX P**

Couple of pads for Hernix, men

**M24643A  
M24643C**

### INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

### CHARACTERISTICS

- Hypoallergenic, thin, soft and comfortable cotton fabric structure
- Reinforced and shaped inner structure
- Under-thigh inner pullers, with calibrated traction
- Closed and shaped front, lined with soft cotton jersey
- Fitted for hernia pads (M24699M)

## Hernia boxer briefs



Sizes:  
**1-8**



Colour:  
**Anthracite grey**



Height:  
**Medium**



Containment:  
**Medium-strong**



Type:  
**Boxer**



Composition:  
**68% CO • 17% PA  
13% EA • 2% PL**



## Female hernia brief



Sizes:  
**1-8**



Colour:  
**Nude**



Height:  
**High**




Containment:  
**Strong**



Type:  
**Female brief**



Composition:  
**82% CO • 13% EA  
5% PA**



**INCLUDED:**  
**M24699M HERNIX P**

Couple of pads for Hernix, men


**M24649**

### INDICATIONS

- Containment for starting hernias and reducible inguinal hernias
- Pre-op (with the insertion of hernia compression pads)
- Post-op
- Recurrence prevention

### CHARACTERISTICS

- Hypoallergenic, thin, soft and comfortable cotton fabric structure adapted for boxer
- Reinforced and shaped inner structure
- Under-thigh inner pullers, with calibrated traction
- Closed and shaped front, lined with soft cotton jersey
- Fitted for hernia pads (M24699M)



**INCLUDED:**  
**M24699W HERNIX P**

Couple of pads for Hernix, female

**M24690**

### INDICATIONS

- Inguinal hernias, ptosis, and prolapse containment
- Post-op
- Hernia recurrence prevention

### CHARACTERISTICS

- Very supportive, hypoallergenic, soft and comfortable cotton fabric structure.
- Reinforced and shaped internal structure for correct support of the abdominal wall
- Soft finishing on the leg area for maximum comfort
- Fitted for hernia pads (M24699W)
- Gusset opening

# CONTAINMENT BRIEFS

hernix 45

## What is their purpose

**Hernix 45:** provide support and containment to the abdominal area following surgery or where there is tissue relaxation.

**Hernix 46 and Hernix 47:** contain and support the abdominal and genital area during sporting activity or in any situation that causes severe strain.

## How do they feel ?

They contain efficiently, are comfortable and breathable. You put them on and forget you are wearing them.

## How do they make you feel?

They are comfortable and breathable while containing efficiently. Easy to wear during the day they will also keep the wearer well protected during sports activities.

## CONTRAINDICATIONS

Not to be used if the hernia is incarcerated, strangulated, irreducible or at an advanced stage of inflammation.

Code	M24645 • M24646 • M24647							
Size	1	2	3	4	5	6	7	8
Hip circumf. cm	83-87	88-92	93-97	98-103	104-109	110-115	116-122	123-130

## Post op containment briefs



Size:  
**1-8**



Colour:  
**Anthracite grey**



Height:  
**Medium**



Containment:  
**Medium-strong**



Type:  
**Brief**



Composition:  
**82% CO • 10% EA  
6% PA • 2% PL**



M24645

## INDICATIONS

- Post-abdominal surgery or in case of tissue relaxation for abdominal wall containment
- Genital support

## CHARACTERISTICS

- Soft, thin, supportive cotton fabric structure
- Reinforced and shaped inner structure for effective containment of the abdominal wall and support of the genital region
- Closed, soft, highly stretchable waistband

## Sports containment briefs



Size:  
**1-8**



Colour:  
**Anthracite grey**



Height:  
**Medium**



Containment:  
**Light**



Type:  
**Brief**



Composition:  
**82% CO • 10% EA  
6% PA • 2% PL**



M24646

### INDICATIONS

- Where preventive genital support is required during daily activities and sport.

### CHARACTERISTICS

- Medium containment cotton fabric structure, thin and soft
- Closed front, with calibrated containment for correct containment of the genital region

## Jockstrap



Size:  
**1-8**



Colour:  
**Anthracite grey**



Height:  
**Low**



Containment:  
**Light**



Type:  
**Jockstrap**



Composition:  
**65% PA • 18% EA  
15% CO • 2% PL**



M24647

### INDICATIONS

- Where preventive genital support is required during daily activities and sport.

### CHARACTERISTICS

- Tailored cotton fabric front for correct containment of the genital region
- Soft elastic under thigh
- Soft finish crotch



# hernix

## compression underwear



[www.orthoservice.com](http://www.orthoservice.com)



**Headquarter: ORTHOSERVICE AG**  
 Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland  
 Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89  
[info@orthoservice.com](mailto:info@orthoservice.com) · [www.orthoservice.com](http://www.orthoservice.com)



**ORTHOSERVICE DEUTSCHLAND GmbH**  
 Flugstraße 8 · D-76532 Baden-Baden · Deutschland  
 Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13  
[info@orthoservice.de](mailto:info@orthoservice.de) · [www.orthoservice.de](http://www.orthoservice.de)

**ORTHOSERVICE POLSKA Sp. z o.o.**  
 ul. Warszawska 416a · 42-209 Częstochowa · Polska  
 Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785  
[info@orthoservice.pl](mailto:info@orthoservice.pl) · [www.orthoservice.pl](http://www.orthoservice.pl)

**OSUK Ltd**  
 Units 47&50 North East BIC  
 Wearfield  
 Sunderland Enterprise Park East  
 Sunderland SR5 2TA  
 Tel +44 (0) 191 51 66 220  
[info@orthoservice.uk](mailto:info@orthoservice.uk) · [www.orthoservice.uk](http://www.orthoservice.uk)

**RO+TEN S.r.l.**  
 Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia  
 Sede operativa e amministrativa:  
 Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia  
 Tel. 0039 0362 18 51 200  
[info@roplusten.com](mailto:info@roplusten.com) · [www.roplusten.com](http://www.roplusten.com)



**ORTHOSERVICE  
 RO+TEN**

Take care feel better