

2-FIX

Stop the trauma,
free the recovery!



ORTHOSERVICE
RO+TEN

Take care feel better

Targeted stability for one or two fingers

**FUNCTIONAL
DESIGN**

DUAL USE:
one or two fingers
depending on the
pathology



**CUSTOM-MOLDABLE
STRUCTURE:**

the aluminum alloy splint
easily adapts to individual
anatomy

When to use 2-FIX?

Proper immobilization is an essential step in treating finger injuries.

2-FIX helps maintain correct alignment, **limiting unwanted movement** and **promoting optimal recovery** conditions.



PALMAR OR DORSAL APPLICATION

for maximum versatility

SIMPLE, QUICK, AND EFFECTIVE ADJUSTMENT SYSTEM

with Velcro® straps
and return rings

SOFT AND COMFORTABLE

even during
prolonged wear

Precise support,
everyday comfort

2-FIX

Two-finger immobilizer
with moldable rigid splint

Code	REF.25002		
Size	1	2	3
Wrist-third phalanx distance (cm)	13-15	15-19	19-23
Color	Dark grey		

Fits right and left

Indications

- Sprains, dislocations, fractures, or other injuries requiring immobilization of one or two fingers
- Pre-operative, post-operative, or post-traumatic use

Contraindications

- Currently no known



SINGLE PATIENT
MULTIPLE USE



RADIO
LUCENT



LOW
PROFILE



CUSTOMIZABLE
FIT



www.orthoservice.com



Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com



ORTHOERVICE DEUTSCHLAND GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. +49 (0) 7221 40 533 0 · Fax 0049 (0) 7221 40 533 20
info@orthoservice.de · www.orthoservice.de

ORTHOERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa · Polska
Tel. +48 (0) 34 340 13 10 · (NIP): 9492246785
info@orthoservice.pl · www.orthoservice.pl

OSUK Ltd
Units 47&50 North East BIC
Wearfield
Sunderland Enterprise Park East
Sunderland SR5 2TA
Tel +44 (0) 191 51 66 220
info@orthoservice.uk · www.orthoservice.uk

RO+TEN S.r.l.
Sede legale: Via Marco De Marchi 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 0362 18 51 200
info@roplusten.com · www.roplusten.com



**ORTHOERVICE
RO+TEN**

Take care feel better